



Switch to Metrodent Brand

Product Guide

Content

Page 6	Bite Blocks
Page 8	Boxing In Wax
Page 10	Gnathic Waxes
Page 12	Modelling Wax - Light Pink
Page 14	Modelling Wax - No.1
Page 16	Modelling Wax - No.2
Page 18	Modelling Wax - No.4
Page 20	Modelling Wax - Metromorphic
Page 22	Modelling Wax - No.6
Page 24	Modelling Wax - College
Page 26	Ribbon Wax
Page 28	Soft Bite Wax
Page 30	Strip Wax
Page 32	Wax Up Sticks
Page 34	Carving Wax Blocks
Page 36	Stick (Cement) Wax
Page 38	Tracing Sticks
Page 40	Dipping Wax
Page 41	Flange Mouthpacks
Page 42	Waxes for Modelling
Page 46	Metrocryn Denture Base - ST Veined
Page 48	Metrocryn Denture Base - Translucent Veined
Page 50	Metrocryn Denture Base - Standard Pink
Page 52	Metrocryn Denture Base - Translucent Non Veined
Page 54	Metrocryn Denture Base - Natural
Page 56	Metrocryn Denture Base - Ethnic
Page 58	Metrocryn Denture Base - Dark V1 Veined
Page 60	Metrocryn Denture Base - Light Veined
Page 62	Metrocryn Denture Base - Clear
Page 64	Metropair - Denture repair acrylic
Page 66	Metrocryn HI (High Impact) Denture Base - Light Pink Veined
Page 68	Orthomet - Orthodontic Acrylic
Page 70	Metrostone - Type 4 Die Stone
Page 71	PCS - Plaster Coating Solution
Page 72	Metrotray - Light curing special trays
Page 73	Matrix Duo - Laboratory Putty
Page 74	Metro Lab Putty - Silicone modelling material

Content

Page 75	Metroform - Custom impression tray acrylic
Page 76	Shellac - Special tray blanks and base plates
Page 77	Gypsolve - Cleaning Fluid
Page 78	Metroclens - Cleaning Fluid
Page 79	Duplicating Gels, Metrogel/Supergel/Plastergel
Page 82	Metrosil Plus - 2 part silicone duplicating material
Page 84	Steribim - Germicidal pumice compound
Page 85	Metroblast - Glass Beads
Page 86	Metroblast - Alumina
Page 88	Metroclave/Hi-Lo - Cold Cure Vessels
Page 90	Wax Wire for Sprues
Page 92	Wax Grid Retentions - Diagonal
Page 94	Wax Grid Retentions - With Holes
Page 96	Wax Grid Retentions - With Mesh
Page 98	Wax Border Strips - With Retentions
Page 100	Wax Edge Strips - Retentions
Page 100	Wax Edge Strips - Perforated Retentions
Page 100	Wax Edge Strips - Round Holes
Page 102	Wax Profile Sticks - Lower Jaw Bars
Page 102	Wax Profile Sticks - Half Round
Page 104	Wax Patterns/Wax Clasp Profiles for - Bonyhard Type A
Page 104	Wax Patterns/Wax Clasp Profiles for - Bonyhard Type B
Page 104	Wax Patterns/Wax Clasp Profiles for - Premolars and Molars
Page 106	Wax Patterns/Wax Clasp Profiles for - Straight/Curved
Page 108	Wax Profile Sticks - Beading Wires
Page 109	Stippled Casting Wax
Page 110	Smooth Casting Wax
Page 111	Preparation Wax
Page 112	Metroglue

Bite Blocks
Boxing In Wax
Gnathic Wax
Modelling Wax
Ribbon Wax
Soft Bite Wax
Strip Wax
Wax Up Sticks

Tried and Trusted



Metrodent Waxes

With over 80 years manufacturing experience Metrodent produce a comprehensive range of dental waxes tried and trusted by dentists and dental technicians both in the UK and around the world.

Using the finest raw materials Metrodent waxes are blended and produced to strict quality control standards giving you a first class British made product.



Made by Metrodent





Bite Blocks

For detailed occlusal records

Metrodent Bite Blocks can be easily adapted to shape and carved with a wax knife. Designed to represent the tooth position and are easy to use. Ideal for a wide range of temperature and climate conditions.

Colours	Pink
	Light Pink
Block Size	114 x 13 x 10 mm
Contents	48 Blocks (8 x 6)
Pack Size	585 g
Order Code - Pink	WAXBB48
Order Code - Light Pink	WAXBB48LP





Boxing In Wax

For boxing impressions prior to casting


Suitable for production of upper and lower jaws, for opposing models as well as individual split-cast bases. It surrounds the casting, can be easily adapted to its form and builds the edge for the model after the pouring of gypsum.

Colour	Red
Sheet Size	300 x 40 x 2.5 mm
Contents	18 Sheets
Pack Size	450 g
Order Code	BIW



Made by Metrodent

Modelling Wax



Gnathic

Differentiate your private work at the bite rim stage Available in Sheet Wax and Bite Blocks

Ivory wax sheets and Ivory wax bite blocks for superior occlusal analysis. More accurately represents the dentate area of the rim. Improves patient communication and the intra oral contrast makes the registration process easier.

Colour	Ivory
Sheet Sizes	183 x 87 x 1.75 mm 183 x 87 x 1.50 mm
Pack Size	500 g
Order Code - 1.75 mm	GWAX500
Order Code - 1.50 mm	GWAX500EXPORT

Bite Blocks	114 x 13 x 10 mm
Contents	48 Blocks (8 x 6)
Pack Size	585 g
Order Code	GWAXBB48



Made by Metrodent

Modelling Wax

The background of the image is a close-up of a red wax block. The wax has a series of parallel, diagonal ridges or striations running from the top-left towards the bottom-right. The color is a vibrant red, and the lighting creates a slight gradient, with the top-right corner appearing slightly lighter than the bottom-left.

Light Pink

A very popular shade

The demand for Metrodent to produce a modelling wax in a light pink shade was overwhelming. Rapidly becoming a best seller our Light Pink modelling wax is based on our No.2 formula. A blend of paraffin and micro-crystalline waxes gives familiar ease of handling and consistent qualities aiding smooth work processing in the laboratory. This desirable colour option offers excellent soft tissue representation at try-in stage and is well received by both patients and clinicians.

Colour	Light Pink
Sheet Sizes	183 x 87 x 1.75 mm 183 x 87 x 1.50 mm
Pack Sizes	500 g 2.5 Kg 20 Kg
Order Code - 1.75 mm x 500 g	LPWAX500
Order Code - 1.75 mm x 2.5 Kg	LPWAX2.5
Order Code - 1.75 mm x 20 Kg	LPWAX20
Order Code - 1.50 mm x 500 g	LPWAX500EXPORT
Order Code - 1.50 mm x 2.5 Kg	LPWAX2.5EXPORT
Order Code - 1.50 mm x 20 Kg	LPWAX20EXPORT



Made by Metrodent



Modelling Wax

CE

No.1

Highest melting point, ideal for carving

Our highest melting point wax originally designed for warmer climates. Utilising a higher melting point paraffin wax giving a longer molecular chain, it is more stable at higher ambient temperatures and therefore gives excellent carving properties. Ideally suited for high end prosthetics where the try-in is extensively contoured and detailed and where the try-in may be in the mouth for an extended period. Perfect for a hot laboratory in the summer months No.1 will be more brittle in lower temperatures.

Colour	Very Pale Pink
Sheet Sizes	183 x 87 x 1.75 mm 183 x 87 x 1.50 mm
Pack Sizes	500 g 2.5 Kg 20 Kg
Order Code - 1.75 mm x 500 g	1WAX500
Order Code - 1.75 mm x 2.5 Kg	1WAX2.5
Order Code - 1.75 mm x 20 Kg	1WAX20
Order Code - 1.50 mm x 500 g	1WAX500EXPORT
Order Code - 1.50 mm x 2.5 Kg	1WAX2.5EXPORT
Order Code - 1.50 mm x 20 Kg	1WAX20EXPORT



Made by Metrodent

Modelling Wax

A close-up photograph of a stack of red modelling wax blocks. The top block is in sharp focus, showing its rectangular shape and slightly textured surface. Below it, several other blocks are visible, but they are increasingly out of focus, creating a sense of depth. The lighting is soft, highlighting the edges and corners of the wax blocks.

No.2

Our most popular blend

This is our most popular modelling wax both in the UK and abroad. It is a blend of paraffin and micro-crystalline waxes and is designed to give good overall handling properties for everyday laboratory use. Pale Pink in colour No.2 is the longest manufactured blend in our range. Ideal for high workflow prosthetic laboratories where its consistency aids smooth production.

Colour

Pale 'Bluish' Pink

Sheet Sizes

183 x 87 x 1.75 mm

183 x 87 x 1.50 mm

183 x 87 x 3.50 mm

Pack Sizes

500 g

2.5 Kg

20 Kg

Order Code - 1.75 mm x 500 g

2WAX500

Order Code - 1.75 mm x 2.5 Kg

2WAX2.5

Order Code - 1.75 mm x 20 Kg

2WAX20

Order Code - 1.50 mm x 500 g

2WAX500EXPORT

Order Code - 1.50 mm x 2.5 Kg

2WAX2.5EXPORT

Order Code - 1.50 mm x 20 Kg

2WAX20EXPORT

Order Code - 3.50 mm x 500 g

DTW

Order Code - 3.50 mm x 2.5 Kg

DTW2.5



Made by Metrodent

Modelling Wax



No.4

Darker Pink variant on our most popular blend

This modelling wax is a colour option of our most popular wax. Its handling properties and indications for use are the same as No.2.

Colour	Darker Pink
Sheet Sizes	183 x 87 x 1.75 mm 183 x 87 x 1.50 mm
Pack Sizes	500 g 2.5 Kg 20 Kg
Order Code - 1.75 mm x 500 g	4WAX500
Order Code - 1.75 mm x 2.5 Kg	4WAX2.5
Order Code - 1.75 mm x 20 Kg	4WAX20
Order Code - 1.50 mm x 500 g	4WAX500EXPORT
Order Code - 1.50 mm x 2.5 Kg	4WAX2.5EXPORT
Order Code - 1.50 mm x 20 Kg	4WAX20EXPORT



Made by Metrodent

Modelling Wax



Metromorphic

A pleasant aroma and easy to work

Metromorphic is a blend of paraffin wax, microcrystalline waxes and natural beeswax which with the omission of resins gives this wax a pleasant aroma. A useful everyday modelling wax.

Colour	Pale 'Bluish' Pink
Sheet Sizes	183 x 87 x 1.75 mm 183 x 87 x 1.50 mm
Pack Sizes	500 g 2.5 Kg 20 Kg
Order Code - 1.75 mm x 500 g	7WAX500
Order Code - 1.75 mm x 2.5 Kg	7WAX2.5
Order Code - 1.75 mm x 20 Kg	7WAX20
Order Code - 1.50 mm x 500 g	7WAX500EXPORT
Order Code - 1.50 mm x 2.5 Kg	7WAX2.5EXPORT
Order Code - 1.50 mm x 20 Kg	7WAX20EXPORT



Made by Metrodent

Modelling Wax

The background of the slide is an abstract composition of soft, wavy lines in various shades of pink and red. In the upper left quadrant, there is a square area with a fine, woven texture, possibly representing a piece of fabric or a specific material. The overall aesthetic is modern and artistic.

No.6

Easy to work with enhanced carving properties

Designed to be easy to work with, No.6 is a blend of paraffin waxes, micro-crystalline waxes, natural beeswax and resins. The addition of resins enhances the carving properties of the wax making it suitable for contouring and characterisation of the denture base from try-in through to fit. Suitable for prosthetic laboratories with high workflow who may need to characterise their try-ins and require an all round reliable modelling wax.

Colour	Pale Pink
Sheet Sizes	183 x 87 x 1.75 mm 183 x 87 x 1.50 mm
Pack Sizes	500 g 2.5 Kg 20 Kg
Order Code - 1.75 mm x 500 g	6WAX500
Order Code - 1.75 mm x 2.5 Kg	6WAX2.5
Order Code - 1.75 mm x 20 Kg	6WAX20
Order Code - 1.50 mm x 500 g	6WAX500EXPORT
Order Code - 1.50 mm x 2.5 Kg	6WAX2.5EXPORT
Order Code - 1.50 mm x 20 Kg	6WAX20EXPORT



Made by Metrodent

Modelling Wax



College

For training and education purposes

Ideal for Dental Schools, Colleges and Universities this wax is made of secondary blends of our primary waxes and varies in colour and handling properties.

Colour	Pink
Sheet Sizes	183 x 87 x 1.75 mm 183 x 87 x 1.50 mm
Pack Sizes	500 g 20 Kg
Order Code - 1.75 mm/500 g	CWAX500
Order Code - 1.75 mm/20 Kg	CWAX20
Order Code - 1.50 mm/500 g	CWAX500EXPORT
Order Code - 1.50 mm/20 Kg	CWAX20EXPORT



Made by Metrodent




Ribbon Wax

For dorsal damming and other intra oral applications

Metrodent Ribbon Wax also known as Periphery Wax is used by dentists and technicians to extend the peripheral edge on an impression tray. High degree of plasticity and exact dorsal dam.

Colour	Red
Contents	126 Strips
Strip Size	210 x 8.0 x 1.50 mm
Pack Size	400g
Order Code	RWAX





Soft Bite Wax

Economical bite registration wax

Metrodent Soft Bite Wax is a proven bite registration material, formulated to provide dentists with a dense yet easily softened bite wax that is stable at room temperature. The thickness has been specially chosen to provide the most accurate registration.

Colour	Red
Contents	24 Strips (6 x 6 sheets)
Strip Size	230 x 12 x 7.5 mm
Pack Size	400g
Order Code	SRBW



Made by Metrodent



Strip Wax

An alternative bite registration wax

Metrodent Strip Wax is a proven bite registration material, formulated to provide dentists with a dense yet easily softened bite wax that is stable at room temperature.

The thickness has been specially chosen to provide the most accurate registration.

Colour	Red
Contents	35 Strips (7 x 5 sheets)
Strip Size	230 x 12 x 7.5 mm
Pack Size	200g
Order Code	SWAX



Made by Metrodent



Wax Up Sticks

For setting up and waxing try-in dentures

Metrodent Wax Up Sticks are specially formulated for the waxing up process of either the bite block or try in stages. It carves easily and boils out cleanly.

Colours

Pink

Contents

Light Pink

Strip Size

64 Sticks

Pack Size

ø 9 x 116 mm

Order Code - Pink

400g

Order Code - Light Pink

WUS

LPWUS



Made by Metrodent



Carving Wax Blocks

Used as a training aid by students

For the study of morphology.

Order Code

Description

Ref: DWS511

Carving Wax Blocks - 30 pieces





little
nm, 12 Sticks per box
W70G

Sticky (Cement) Wax

For plaster, metals and resins

Ideal for sticking to plaster, metals and resin

Order Code

Description

Ref: SW70G

Sticky (Cement) Wax - 70g

Ref: SW500G

Sticky (Cement) Wax - 50g





Tracing Sticks

Impression compound

For adjusting the periphery of impression trays

Order Code

Description

Ref: COM621

Tracing Sticks - 15 per box





Dipping Wax

Granular

Cost effective way of producing crown and bridge copings

Order Code	Description
Ref: MDWG1KG	Dipping Wax - 1 Kg



Flange Mouthpacks

Non Sterile

Absorbent sponges to pack the airway in order to prevent inhalation/ ingestion of debris, etc.

Autoclave before use.

Order Code	Description
Ref: MPACK	Flange Mouth Packs - 12 per pack



Cervical Wax

Crown Wax

Undercut Wax

For modelling

Order Code

Description

Ref: W24530

Cervical Wax - Red 45g

Ref: W24511

Crown Wax - Blue (Extra Hard) 45g

Ref: W23002

Undercut Wax (Block Out) - Pink 50g





Metrocyl Universal
Metrocyl Rapid Cure
Metrocyl HI
Metropair
Orthomet
Metroform

Tried and Trusted



Metrodent Denture Base Acrylics

Over decades Metrodent have developed a range of denture base acrylics for both our UK and Export markets. Metrocryl Rapid Cure Denture Base is available in 9 shade options including clear and ethnic shades. In addition we have an orthodontic resin, Orthomet, for the “spray” technique available in clear and tinted pink.

**All our heat cure resins conform to Standard
BS EN ISO20795-1: 2008**




Made by Metrodent

Metrocryn Denture Base



ST Veined



Pink Translucent Veined

Our most popular acrylic formulated in response to market demand to match the leading brand product. A silky dough is produced with exceptionally good flow properties resulting in reduced flash thickness.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code	Description
Ref: U1SV	Metrocyl - Universal ST Veined 1 Kg
Ref: U3SV	Metrocyl - Universal ST Veined 3 Kg
Ref: U12SV	Metrocyl - Universal ST Veined 12 Kg
Ref: R1SV	Metrocyl - Rapid Cure ST Veined 1 Kg
Ref: R3SV	Metrocyl - Rapid Cure ST Veined 3 Kg
Ref: R12SV	Metrocyl - Rapid Cure ST Veined 12 Kg
Ref: LU500	Metrocyl - Universal Liquid 500 ml
Ref: LU1	Metrocyl - Universal Liquid 1 Litre
Ref: MP1SV	Metropair - Universal ST Veined 1 Kg
Ref: MP3SV	Metropair - Universal ST Veined 3 Kg
Ref: MP12SV	Metropair - Universal ST Veined 12 Kg
Ref: MPL500	Metropair - Universal Liquid 500 ml
Ref: MPL1	Metropair - Universal Liquid 1 Litre




Made by Metrodent



Metrocryn Denture Base



Translucent Veined



Medium Pink Veined

A popular translucent denture base with veining. Slightly whiter or more opaque than either our Translucent or Light Veined acrylic.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code	Description
Ref: U1TV	Metrocryn - Universal Translucent Veined 1 Kg
Ref: U3TV	Metrocryn - Universal Translucent Veined 3 Kg
Ref: U12TV	Metrocryn - Universal Translucent Veined 12 Kg
Ref: R1TV	Metrocryn - Rapid Cure Translucent Veined 1 Kg
Ref: R3TV	Metrocryn - Rapid Cure Translucent Veined 3 Kg
Ref: R12TV	Metrocryn - Rapid Cure Translucent Veined 12 Kg
Ref: LU500	Metrocryn - Universal Liquid 500 ml
Ref: LU1	Metrocryn - Universal Liquid 1 Litre
Ref: MP1TV	Metropair - Translucent Veined 1 Kg
Ref: MP3TV	Metropair - Translucent Veined 3 Kg
Ref: MP12TV	Metropair - Translucent Veined 12 Kg
Ref: MPL500	Metropair - Liquid 500 ml
Ref: MPL1	Metropair - Liquid 1 Litre




Made by Metrodent



Metrocryn Denture Base



Standard Pink



Natural Pink

With greater opacity this acrylic has no veining. A lustrous natural pink mimicking healthy gingival mucosa, un-veined.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code	Description
Ref: U1S	Metrocyl - Universal Standard Pink 1 Kg
Ref: U3S	Metrocyl - Universal Standard Pink 3 Kg
Ref: U12S	Metrocyl - Universal Standard Pink 12 Kg
Ref: R1S	Metrocyl - Rapid Cure Standard Pink 1 Kg
Ref: R3S	Metrocyl - Rapid Cure Standard Pink 3 Kg
Ref: R12S	Metrocyl - Rapid Cure Standard Pink 12 Kg
Ref: LU500	Metrocyl - Universal Liquid 500 ml
Ref: LU1	Metrocyl - Universal Liquid 1 Litre
Ref: MP1S	Metopair - Standard Pink 1 Kg
Ref: MP3S	Metopair - Standard Pink 3 Kg
Ref: MP12S	Metopair - Standard Pink 12 Kg
Ref: MPL500	Metopair - Liquid 500 ml
Ref: MPL1	Metopair - Liquid 1 Litre




Made by Metrodent



Metrocryn Denture Base



Translucent Non Veined



Light Gingival Pink

Translucent denture base, translucency allows the patient's own gingival tones to "colour" the denture base.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code

Description

Ref: U1T	Metrocyl - Universal Translucent Non Veined 1 Kg
Ref: U3T	Metrocyl - Universal Translucent Non Veined 3 Kg
Ref: U12T	Metrocyl - Universal Translucent Non Veined 12 Kg
Ref: R1T	Metrocyl - Rapid Cure Translucent Non Veined 1 Kg
Ref: R3T	Metrocyl - Rapid Cure Translucent Non Veined 3 Kg
Ref: R12T	Metrocyl - Rapid Cure Translucent Non Veined 12 Kg
Ref: LU500	Metrocyl - Universal Liquid 500 ml
Ref: LU1	Metrocyl - Universal Liquid 1 Litre

Ref: MP1T	Metropair - Translucent Non Veined 1 Kg
Ref: MP3T	Metropair - Translucent Non Veined 3 Kg
Ref: MP12T	Metropair - Translucent Non Veined 12 Kg
Ref: MPL500	Metropair - Liquid 500 ml
Ref: MPL1	Metropair - Liquid 1 Litre



Made by Metrodent



Metrocryn Denture Base



Natural

Medium Pink

A pink opaque denture base mimicking certain gingival tissues. Medium pink with a very subtle orange tint and good opacity.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code

Description

Ref: U1N	Metrocyl - Universal Natural 1 Kg
Ref: U3N	Metrocyl - Universal Natural 3 Kg
Ref: U12N	Metrocyl - Universal Natural 12 Kg
Ref: R1N	Metrocyl - Rapid Cure Natural 1 Kg
Ref: R3N	Metrocyl - Rapid Cure Natural 3 Kg
Ref: R12N	Metrocyl - Rapid Cure Natural 12 Kg
Ref: LU500	Metrocyl - Universal Liquid 500 ml
Ref: LU1	Metrocyl - Universal Liquid 1 Litre

Ref: MP1N	Metropair - Natural 1 Kg
Ref: MP3N	Metropair - Natural 3 Kg
Ref: MP12N	Metropair - Natural 12 Kg
Ref: MPL500	Metropair - Liquid 500 ml
Ref: MPL1	Metropair - Liquid 1 Litre




Made by Metrodent



Metrocryn Denture Base



Ethnic



Dark Purple shade of Pink veined

Developed to match the gingival tissues of patients with afro-caribbean heritage.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code

Description

Ref: U1E Metrocyl - Universal Ethnic 1 Kg
Ref: U3E Metrocyl - Universal Ethnic 3 Kg
Ref: R1E Metrocyl - Rapid Cure Ethnic 1 Kg
Ref: R3E Metrocyl - Rapid Cure Ethnic 3 Kg
Ref: LU500 Metrocyl - Universal Liquid 500 ml
Ref: LU1 Metrocyl - Universal Liquid 1 Litre

Ref: MP1E Metropair - Ethnic 1 Kg
Ref: MP3E Metropair - Ethnic 3 Kg
Ref: MPL500 Metropair - Liquid 500 ml
Ref: MPL1 Metropair - Liquid 1 Litre




Made by Metrodent



Metrocryn Denture Base





Dark V1 Veined

A more reddish Pink veined

Developed originally for export to our East Asian markets and offering a variation in gingival tissue matching.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code

Description

Ref: U1DV	Metrocryn - Universal Dark V1 Veined 1 Kg
Ref: U3DV	Metrocryn - Universal Dark V1 Veined 3 Kg
Ref: R1DV	Metrocryn - Rapid Cure Dark V1 Veined 1 Kg
Ref: R3DV	Metrocryn - Rapid Cure Dark V1 Veined 3 Kg
Ref: LU500	Metrocryn - Universal Liquid 500 ml
Ref: LU1	Metrocryn - Universal Liquid 1 Litre
Ref: MP1DV	Metropair - Dark V1 Veined 1 Kg
Ref: MP3DV	Metropair - Dark V1 Veined 3 Kg
Ref: MPL500	Metropair - Liquid 500 ml
Ref: MPL1	Metropair - Liquid 1 Litre




Made by Metrodent



Metrocryn Denture Base





Light Veined

Light Translucent Pink Veined

Similar to our most popular ST Veined acrylic but in a lighter more subtle pink suitable for patients with paler gingival tissues.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Order Code

Description

Ref: U1LV	Metrocyl - Universal Light Veined 1 Kg
Ref: U3LV	Metrocyl - Universal Light Veined 3 Kg
Ref: U12LV	Metrocyl - Universal Light Veined 12 Kg
Ref: R1LV	Metrocyl - Rapid Cure Light Veined 1 Kg
Ref: R3LV	Metrocyl - Rapid Cure Light Veined 3 Kg
Ref: R12LV	Metrocyl - Rapid Cure Light Veined 12 Kg
Ref: LU500	Metrocyl - Universal Liquid 500 ml
Ref: LU1	Metrocyl - Universal Liquid 1 Litre
Ref: MP1LV	Metropair - Translucent Non Veined 1 Kg
Ref: MP3LV	Metropair - Translucent Non Veined 3 Kg
Ref: MP12LV	Metropair - Translucent Non Veined 12 Kg
Ref: MPL500	Metropair - Liquid 500 ml
Ref: MPL1	Metropair - Liquid 1 Litre




Made by Metrodent



Metrocryn Denture Base



Clear



Clear

Ideal for clear palates, TMJ appliances, orthodontic appliances, surgical guides etc. Using rapid cure gives an exceptionally clear acrylic whereas universal can sometimes have a slight opacity.

Universal

- ✓ Snap stage 15 minutes
- ✓ Work time 60 minutes at 23°C
- ✓ Can be quick cured or processed overnight

Rapid Cure

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes
- ✓ Use of rapid cure gives a very clear acrylic*

Order Code	Description
Ref: U1C	Metrocryn - Universal Clear 1 Kg
Ref: U3C	Metrocryn - Universal Clear 3 Kg
Ref: U12C	Metrocryn - Universal Clear 12 Kg
Ref: R1C	Metrocryn - Rapid Cure Clear 1 Kg
Ref: R3C	Metrocryn - Rapid Cure Clear 3 Kg
Ref: R12C	Metrocryn - Rapid Cure Clear 12 Kg
Ref: LC500	Metrocryn - Universal Clear Liquid 500 ml
Ref: LC1	Metrocryn - Universal Clear Liquid 1 Litre
Ref: MP1C	Metropair - Clear 1 Kg
Ref: MP3C	Metropair - Clear 3 Kg
Ref: MP12C	Metropair - Clear 12 Kg
Ref: MPL500	Metropair - Liquid 500 ml
Ref: MPL1	Metropair - Liquid 1 Litre

*The use of Universal liquid with clear acrylic will result in yellowing of the cured acrylic and we recommend always using our specially formulated clear liquid.



Made by Metrodent





Metropair

Denture Repair Acrylic

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**

Order Code	Description
Ref: MP1SV	Metropair - Universal ST Veined 1 Kg
Ref: MP3SV	Metropair - Universal ST Veined 3 Kg
Ref: MP12SV	Metropair - Universal ST Veined 12 Kg
Ref: MP1TV	Metropair - Translucent Veined 1 Kg
Ref: MP3TV	Metropair - Translucent Veined 3 Kg
Ref: MP12TV	Metropair - Translucent Veined 12 Kg
Ref: MP1S	Metropair - Standard Pink 1 Kg
Ref: MP3S	Metropair - Standard Pink 3 Kg
Ref: MP12S	Metropair - Standard Pink 12 Kg
Ref: MP1T	Metropair - Translucent Non Veined 1 Kg
Ref: MP3T	Metropair - Translucent Non Veined 3 Kg
Ref: MP12T	Metropair - Translucent Non Veined 12 Kg
Ref: MP1N	Metropair - Natural 1 Kg
Ref: MP3N	Metropair - Natural 3 Kg
Ref: MP12N	Metropair - Natural 12 Kg
Ref: MP1E	Metropair - Ethnic 1 Kg
Ref: MP3E	Metropair - Ethnic 3 Kg
Ref: MP1DV	Metropair - Dark V1 Veined 1 Kg
Ref: MP3DV	Metropair - Dark V1 Veined 3 Kg
Ref: MP1LV	Metropair - Translucent Non Veined 1 Kg
Ref: MP3LV	Metropair - Translucent Non Veined 3 Kg
Ref: MP12LV	Metropair - Translucent Non Veined 12 Kg
Ref: MP1C	Metropair - Clear 1 Kg
Ref: MP3C	Metropair - Clear 3 Kg
Ref: MP12C	Metropair - Clear 12 Kg
Ref: MPL500	Metropair - Liquid 500 ml
Ref: MPL1	Metropair - Liquid 1 Litre



Made by Metrodent



Metrocryn Denture Base



HI (High Impact)

Light Pink Veined

Metrocyl HI, when cured, forms a virtually unbreakable prosthesis. The impact strength of notched and surface defected specimens (according to ISO 179) gives values of 2.55 Kj/m² and 9.55 Kj/m² respectively.

Polymer

- ✓ Snap stage 15 minutes
- ✓ Work time 20 minutes at 23°C

Monomer

- ✓ Only use with Metrocyl HI

Order Code

Description

Ref: HI1	Metrocyl HI - Polymer 1 Kg
Ref: HI3	Metrocyl HI - Polymer 3 Kg
Ref: HI12	Metrocyl HI - Polymer 12 Kg
Ref: HIL500	Metrocyl HI - Monomer 500 ml
Ref: HIL1	Metrocyl HI - Monomer 1 Litre

Impact strength is different from tensile or flexural strength and will prevent breakage of the denture base due to impact shocks for example if the denture is dropped on a hard surface like a porcelain sink. High impact acrylic is not a substitute for a CoCr denture base where a patient's bite or biting forces would indicate this prescription. In thin cross sections, such as the replacement of a single tooth such as a central or lateral incisor in a close bite situation, high impact acrylics are also no more suitable than a conventional PMMA denture base. Repairs – Should the need arise to repair a denture fabricated from high impact denture base it should be re-processed using heat cure high impact resin to maintain its high impact properties. The use of a standard "cold cure" resin will mean the denture is no longer considered to have greater impact resistance. Unlike heat cure denture base, the high impact powder does not absorb monomer very well due to it being rubber modified. To this end the powder should be weighed and the volume of liquid measured according to the instructions. The powder will need to be worked into the liquid using a metal spatula to ensure the correct amount of powder is used in the mix. Use of the 'pepper and salt' technique is likely to give incorrect polymerisation.




Made by Metrodent



Orthodontic Acrylic



Orthomet



Clear or Tinted Pink

Orthomet is used in the “spray” technique for the manufacture of Orthodontic appliances. With the use of a pressure flask such as the Metroclave or Hi-Lo, the material will be cured in approximately 10 minutes showing no porosity or shrinkage.

Order Code	Description
Ref: OP500	Orthomet - Polymer 500 g
Ref: OP1	Orthomet - Polymer 1 Kg
Ref: OP3	Orthomet - Polymer 3 Kg
Ref: OP10	Orthomet - Polymer 10 Kg
Ref: OLC250	Orthomet - Clear Monomer 250 ml
Ref: OLC500	Orthomet - Clear Monomer 500 ml
Ref: OLC1	Orthomet - Clear Monomer 1 Litre
Ref: OLT250	Orthomet - Tinted Pink Monomer 250 ml
Ref: OLT500	Orthomet - Tinted Pink Monomer 500 ml
Ref: OLT1	Orthomet - Tinted Pink Monomer 1 Litre



Made by Metrodent






Metrostone

Type 4 Die Stone

A fine grain structure ensures a very smooth creamy mix. Exceptional strength and durability. Setting time between 8 and 12 minutes.

Order Code	Description
MSTONEB	Metrostone - Blue 22.68 Kg
MSTONEI	Metrostone - Ivory 22.68 Kg
MSTONEW	Metrostone - White 22.68 Kg
MSTONEP	Metrostone - Pink 22.68 Kg
MSTONEG	Metrostone - Green 22.68 Kg
MSTONEPEACH	Metrostone - Peach 22.68 Kg



Plaster Coating Solution

Cold Mould Seal

Premium grade formaldehyde free alginate solution for insulation of denture base resins from gypsum moulds and casts when making dentures, acrylic crowns etc.

Order Code	Description
Ref: PCS1	PCS (Plaster Coating Solution) - 1 Litre
Ref: PCS2	PCS (Plaster Coating Solution) - 2 Litres
Ref: PCS5	PCS (Plaster Coating Solution) - 5 Litres
Ref: PCS25	PCS (Plaster Coating Solution) - 25 Litres
Ref: PCS200	PCS (Plaster Coating Solution) - 200 Litres



Metrotray

Light curing special tray material

Light curing hybrid composite resin for making individual tray blanks.

Order Code	Description
Ref: DP08-00900	Metrotray Pink - 50 pack
Ref: DP08-00901	Metrotray Blue - 50 pack



Matrix Duo

Laboratory Putty

Addition cured 2 part silicone lab putty.
Sh Hardness 8 minutes (75/96), 24 hours (82/96)

Order Code	Description
Ref: MP2	Matrix Duo - Laboratory Putty 2 x 1.5 Litre



Metro Lab Putty

Silicone modelling material

Supplied as a base material and catalyst it is easily mixed by hand. Lab putty is routinely used for the production of copy denture moulds and denture repairs. Increasingly lab putty is being used for processing denture repairs. After sticking the parts of the broken denture together the repair model can be moulded quickly into the denture and placed on the bench to set.

The fit surface is accurately reproduced before and after processing. The use of cold mould seal is not required. If for any reason the repair cannot be processed immediately, there is no danger of bubbles appearing in the cold cure.

When repairing a line crack it is advisable to make a complete model rather than one just over the crack as the denture can sometimes fracture when grinding out the repair. Lab putty is also used for matrices in the construction of partial dentures, for positioning backings, cold curing CoCr saddles etc.

Order Code	Description
Ref: MLP	Metro Lab Putty - 5 Kg
Ref: MLPC	Metro Lab Putty - Activator 40 ml
3 to 4 activators recommended for 5 Kg of Metro Lab Putty	



Metroform

Custom impression tray acrylic

A high quality self cure material with excellent working properties and economic in use.

Blue - standard time.
Dough time 2 to 3 minutes.

Green - long work time.
Dough time 4 to 5 minutes.

Order Code	Description
Ref: MTFB1	Metroform Blue - Polymer 1 Kg
Ref: MTFB3	Metroform Blue - Polymer 3 Kg
Ref: MTFB12	Metroform Blue - Polymer 12 Kg
Ref: MTFG1	Metroform Green - Polymer 1 Kg
Ref: MTFG3	Metroform Green - Polymer 3 Kg
Ref: MTFG12	Metroform Green - Polymer 12 Kg
Ref: MTFL500	Metroform - Monomer 500 ml
Ref: MTFL1	Metroform - Monomer 1 Litre



Shellac

Special Tray Blanks and Base Plates

Special Tray blanks are a firm favourite with Dentists and Technicians. Easily adapted to produce individual custom trays for the most accurate impressions.

Use Shellac Base plates when a very rigid base is needed to ensure maximum accuracy at the try-in stage of a denture.

Order Code	Description
Ref: PAT112	Shellac Special Tray Blanks - Green, XL Plain Upper (2.25 mm deep) x 72
Ref: PAT122	Shellac Special Tray Blanks - Green, XL Perforated (2.25 mm deep) x 72
Ref: BPL003/P	Shellac (Perfecta) Base Plates - Pink, Upper (1.25 mm deep) x 100



Gypsolve

Cleaning Fluid

Gypsolve is a clear non toxic, non corrosive liquid which can safely be used for the rapid removal of residual plaster. Also suitable for use in ultrasonic cleaners.

Order Code

Description

Ref: GG5

Gypsolve - Green 5 Litres

Ref: GG25

Gypsolve - Green 25 Litres

Ref: GC5

Gypsolve - Clear 5 Litres

Ref: GC25

Gypsolve - Clear 25 Litres



Made by Metrodent

Metroclens

Cleaning Fluid

Metroclens is a tray cleaner, plaster remover and hard surface cleaner. Used in varying diluted forms gives Metroclens solution optimum cleaning properties. Can be used in an ultrasonic cleaner to speed up the cleaning process. Must not be used for aluminium or aluminium alloys.

Order Code	Description
Ref: MCLENS1	Metroclens - 1 Litres



Made by Metrodent

Duplicating Gels

Metrogel

Supergel

Plastergel

Metrogel - successful aid in the partial denture technique.

Supergel - for the cobalt chrome partial denture laboratory.

Plastergel - is ideal for the production of quality duplicate models made from plasters and dental stones.

Order Code

Description

Ref: MGEL

Metrogel - Blue 5 Kg

Ref: MGEL/C

Metrogel - Clear 5 Kg

Ref: SGEL

Supergel - Green 5 Kg

Ref: SGEL/C

Supergel - Clear 5 Kg

Ref: PGEL

Plastergel - Red 5 Kg

Ref: PGEL/C

Plastergel - Clear 5 Kg



Made by Metrodent





$A + B =$



Metrosil Plus

2 part silicone duplicating material

Metrosil is a two part addition cured silicone for the accurate duplication of master models. Mixed in a 1:1 ratio the material is blended to a uniform colour before pouring. Compared to reversible gels, silicone duplication has always been considered expensive. However, the growth of private and implant work has placed demands on the dental laboratory for high levels of accuracy and precision. Metrosil has excellent dimensional stability at an economical price.

Order Code	Description
Ref: MSDPLUS	Metrosil Plus - 1 Kg Base and 1 Kg Catalyst





Steribim

Germicidal pumice compound

Destroys all germs and the germicidal effect lasts down to the last gram. Special additives make it kind to the skin and prevent the pumice from going sour in the polishing trough. Dust free.

Order Code

Description

Ref: STB4

Steribim Germicidal Pumice - 4 Kg

Ref: STB12

Steribim Germicidal Pumice - 12 Kg



Made by Metrodent

Metroblast Glass Beads

Blasting Material

In the manufacture of the glass beads alkali-lime glass is broken and prescreened. Subsequently, the glass fragments are formed in a furnace in the heat flow to balls. Glass beads are a mineral abrasive.

Advantages

- ✓ Glass beads for blast polishing
- ✓ No surface metal loss

Order Code	Description
Ref: MBGB50/10	Metroblast Glass Beads 50µm - 10 Kg
Ref: MBGB50/25	Metroblast Glass Beads 50µm - 25 Kg



Metroblast Alumina

Blasting Material

Ultra pure, Iron Free Aluminium Oxide for blast polishing with no surface metal loss. Metrodent blasting material is a white fused alumina manufactured by melting a special grade of calcined alumina (aluminium oxide) in an electric arc furnace using the block method; it is ultra pure aluminium oxide. An iron-free fine abrasive/refractory material for the toughest demands.

Advantages

- ✓ Ultra pure, Iron free aluminium oxide
- ✓ No surface metal loss

Order Code	Description
Ref: MBAO250/10	Metroblast Alumina 250µm - 10 Kg
Ref: MBAO110/10	Metroblast Alumina 150µm - 10 Kg
Ref: MBAO50/10	Metroblast Alumina 50µm - 10 Kg
Ref: MBAO250/25	Metroblast Alumina 250µm - 25 Kg
Ref: MBAO110/25	Metroblast Alumina 150µm - 25 Kg
Ref: MBAO50/25	Metroblast Alumina 50µm - 25 Kg



Cold Cure Vessels



Metroclave and Hi-Lo

The simple and reliable solution

Metrodent have been making the Metroclave and Hi-Lo Cold Cure Vessels for over 30 years. We think they are the best, and many Technicians agree with us. Part of their success is based on their simplicity. To process cold cure resins two things are required heat and pressure. Using the Metroclave and Hi-Lo the heat is supplied by filling the vessels with warm water and the pressure is applied by screwing down the lid and indicated on the gauge. That's it!

Advantages

- ✓ No compressor required
- ✓ No heating plate needed
- ✓ Does not tie-up the flask press whilst in use
- ✓ Requires little maintenance
- ✓ Gives years of service
- ✓ Inexpensive

Metroclave: The lid is recessed in a well in the body of the vessel, this allows it to be used on a bench away from the sink or drainer. If the Metroclave is filled to just the right level, the well allows the lid to be fitted without spillage. Curing time 6 to 15 minutes according to the water temperature. Holds up to 5 models. Dimensions: ø 160 x 90 mm.

Hi-Lo: This is smaller of the two vessels and has a capacity for up to three models. It is ideal for the smaller laboratory or department or even for occasional use in a dental surgery. Curing time 6 to 15 minutes according to the water temperature. Dimensions; ø 110 x 110 mm.

Order Code	Description
Ref: MTC	Metrotroclave cold cure vessel
Ref: HLA	Hi-Lo cold cure vessel





Wax Wire for Sprues

Easy to shape and burn-out without residue

Order Code

Description

Ref: W12120

ø 2mm Wax Wire for Sprues - Green 2/3 hard, 250g

Ref: W12125

ø 2.5mm Wax Wire for Sprues - Green 2/3 hard, 250g

Ref: W12130

ø 3mm Wax Wire for Sprues - Green 2/3 hard, 250g

Ref: W12135

ø 3.5mm Wax Wire for Sprues - Green 2/3 hard, 250g

Ref: W12140

ø 4mm Wax Wire for Sprues - Green 2/3 hard, 250g

Ref: W12150

ø 5mm Wax Wire for Sprues - Green 2/3 hard, 250g



Wax Grid Retentions



Diagonal

For maxillary partial denture frames

Wax grid retentions permit the simple and effective shaping of retentions to total or partial dentures. They guarantee a high level of security in the connection between the resin and the partial denture plate. The large grid retentions facilitate very economical use of material.

Colour: Red, 10 sheets
Dimensions 150 x 75 mm

Order Code	Description
Ref: W13030	Wax Grid Retention - Diagonal



Wax Grid Retentions



with holes

For upper jaw

Wax grid retentions permit the simple and effective shaping of retentions to total or partial dentures. They guarantee a high level of security in the connection between the resin and the partial denture plate. The large grid retentions facilitate very economical use of material.

Colour: Red, 20 sheets
Dimensions 65 x 75 mm

Order Code

Description

Ref: W13020

Wax Grid Retention - Diagonal



Wax Grid Retentions



with mesh

For maxillary partial denture frames

Wax grid retentions permit the simple and effective shaping of retentions to total or partial dentures. They guarantee a high level of security in the connection between the resin and the partial denture plate. The large grid retentions facilitate very economical use of material.

Small Mesh - Colour: Red, 25 sheets
Dimensions 40 x 60 mm

Large Mesh - Colour: Red, 25 sheets
Dimensions 100 x 100 mm

Order Code

Description

Ref: W13000

Wax Grid Retention - Small Mesh Red 25 Sheets

Ref: W13010

Wax Grid Retention - Large Mesh Red 25 Sheets



Wax border strips



with retentions


Time saving

Time savings when modelling upper-jaw frames with large saddles. A great advantage is that the border strip can easily be shaped as desired since the size can be varied by trimming the tips of the retentions. The wax is so supple that it can be shaped easily and reliably as required onto the duplicate model.

Order Code	Description
Ref: W14010	Wax Border Strips - Small (100 mm) with retentions Red 15 pieces
Ref: W14020	Wax Border Strips - Large (130 mm) with retentions Red 15 pieces



Wax Edge Strips with



retentions

perforated retentions

round holes


For lower jaw partial denture frames

For the secure attachment of plastic saddles to lower partial dentures.
Wax edge strips with retentions - Colour: Red, 25 pieces (120 mm length)
Wax edge strips with perforated retentions - Colour: Red, 25 pieces (125 mm length)
Wax edge strips with round holes - Colour: Red, 25 pieces (125 mm length)

Order Code	Description
Ref: W14110	Wax Edge Strips With Retentions (120 mm) Red 25 pieces
Ref: W14120	Wax Edge Strips With Perforated Retentions (125 mm) Red 25 pieces
Ref: W13040	Wax Edge Strips With Round Holes (135 mm) Red 25 pieces



Wax Profile Sticks



Lower Jaw Bars, Half Round

Extremely versatile

Can be used for secondary spruing to clasps, forming continuous clasps or minor connectors.

Lower Jaw Bars - Colour: Green, (65g) Dimensions 145 mm length

Half Round (continuous clasps) - Colour: Green 50g Dimensions 1.8 x 0.9 x 145 mm, Dimensions 2.0 x 1.0 x 145 mm

Order Code	Description
Ref: W12251	Wax Profile Sticks Lower Jaw Bars - 4.0 x 1.7 mm Green 65g
Ref: W12252	Wax Profile Sticks Lower Jaw Bars - 4.3 x 2.3 mm Green 65g
Ref: W12240	Wax Profile Sticks Lower Jaw Bars - 5.0 x 1.3 mm Green 65g
Ref: W12250	Wax Profile Sticks Lower Jaw Bars - 4.0 x 2.0 mm Green 65g
Ref: W12212	Wax Profile Sticks Half Round (Continuous Clasps) - 1.8 mm Green 50g
Ref: W12213	Wax Profile Sticks Half Round (Continuous Clasps) - 2.0 mm Green 50g



Wax Patterns/ Wax Clasp Profiles for



Bonyhard Type A Bonyhard Type B Premolars and Molars

Extremely versatile

These preshaped clasp profiles simplify modelling and help save time. The shape of the profiles enables a large number of variations by shortening the wax form.

Order Code

Description

Ref: W50020

Wax Patterns/Wax Clasp Profiles for Bonyhard Type A - Green 200 clasps

Ref: W50021

Wax Patterns/Wax Clasp Profiles for Bonyhard Type B - Green 200 clasps

Ref: W50010


Wax Patterns/Wax Clasp Profiles for Premolars - Green 200 clasps

Ref: W50000

Wax Patterns/Wax Clasp Profiles for Molars - Green 200 clasps



Wax Patterns/ Ring Clasp Profiles



Straight and Curved

Standard clasps

Slender standard clasp profiles for molars.
Colour: Green
(10 sheets) 200 clasps

Order Code	Description
Ref: W50031	Wax Patterns/Ring Clasp Profiles - Straight Green 200 clasps
Ref: W50030	Wax Patterns/Ring Clasp Profiles - Curved Green 200 clasps



Wax Profile Sticks



Beading Wires

Neat and accurate finish to saddle areas

Order Code

Description

Ref: W12201

Wax Profile Sticks - Beading Wire, Green 25g \varnothing 0.8 mm x 145 mm

Ref: W12203

Wax Profile Sticks - Beading Wire, Green 25g \varnothing 1.0 mm x 145 mm

Ref: W12207

Wax Profile Sticks - Beading Wire, Green 25g \varnothing 1.5 mm x 145 mm

Ref: W12208

Wax Profile Sticks - Beading Wire, Green 25g \varnothing 2.0 mm x 145 mm



Stippled Casting Wax

For modelling the bases of upper partial dentures

Tried and tested wax for modelling the bases of upper partial dentures. Can be easily adapted and adheres firmly to the investment model with no additional wax adhesive. The stippled casting wax is available in three different surface textures – from fine to coarse – and allows customisation of the surface shape as required by the practitioner. The individual stippling of the cast partial denture base facilitates the gripping of food and reduces the foreign body sensation for the patient's tongue. Colour: Green (15 sheets)
Dimensions 150 x 75 mm

Order Code	Description
Ref: W11035	Stippled Casting Wax - Fine 0.35 mm
Ref: W11040	Stippled Casting Wax - Fine 0.40 mm
Ref: W11050	Stippled Casting Wax - Fine 0.50 mm
Ref: W11060	Stippled Casting Wax - Fine 0.60 mm
Ref: W11340	Stippled Casting Wax - Medium 0.40 mm
Ref: W11140	Stippled Casting Wax - Coarse 0.40 mm
Ref: W11150	Stippled Casting Wax - Coarse 0.50 mm
Ref: W11160	Stippled Casting Wax - Coarse 0.60 mm



Smooth Casting Wax

Simple crease free adaptation

Simple, crease-free adaptation. Adheres firmly to the investment model and burns out leaving no residue. The high transparency of the wax makes for optimal clarity of the construction markings on the master model and saves unnecessary, time-consuming corrections to the wax-up.

Colour: Green (15 sheets)

Dimensions 175 x 80 mm

Order Code	Description
Ref: W10030	Smooth Casting Wax - 0.30 mm
Ref: W10040	Smooth Casting Wax - 0.40 mm
Ref: W10050	Smooth Casting Wax - 0.50 mm
Ref: W10060	Smooth Casting Wax - 0.60 mm
Ref: W10070	Smooth Casting Wax - 0.70 mm



Preparation Wax

For the partial denture technique

The preparation wax is exceptionally malleable, allowing it to be adapted to the master model perfectly and with firm adhesion, which saves having to use an additional wax adhesive. The exemplary shape retention and edge strength of the preparation wax, with a high solidification point of approx. 70°C, mean that it can be used with duplicating hydrocolloid at working temperatures of 55°C. Simple removal from the master model following duplication rounds off the user friendly working characteristics perfectly.

Colour: Red (15 sheets)

Dimensions 175 x 80 mm

Order Code	Description
Ref: W11750	Preparation Wax - Fine 0.50 mm
Ref: W11760	Preparation Wax - Fine 0.60 mm
Ref: W11770	Preparation Wax - Fine 0.70 mm



Metroglue

Cyanoacrylate Adhesive (Super Glue)

General purpose super glue.
Bonds in seconds. 20 g

Spray activator that accelerates the adhesive properties of the glue.

Order Code

Description

Ref: HGLUE

Metroglue Adhesive - 20g

Ref: VA200

Metroglue Accelerator - 200 ml





Metrodent.com

