

# Blasters

## Have you seen the light?

Introducing the new range of BEGO Sandblasters, the first with LED technology



All new BEGO Sandblasters are equipped with advanced LED technology, 100% more light (5000 lux) for shadow free sandblasting, this ensures maximum light in the areas where the dental technician needs it most, both in the blasting chamber and on the handpiece.

Brighter, more robust and extremely cost effective, BEGO's new range of Sandblasters gives you the reliability and confidence to prepare more professionally and at the highest level of technology.

### EasyBlast

only £1,129.12\* Ref: 26385

- Highly integrated LED lighting in the handpiece for an optimal shadow-free view during preparation, even intracoronally
- It is possible to detect cracks in the ceramic during sandblasting with the aid of the LED technology – earliest possible quality management and avoidance of complaints at a later stage
- LED lighting in the sandblasting chamber provides a service life ten times higher than that of conventional lamps for cost-effective sandblasting
- The new EasyBlast can be upgraded to a 4-chamber sandblaster using 2 other sandblasting modules – you save investment in another unit
- Intelligent colour coding system on the tanks and in the sandblasting chamber with the aid of LED lights, making it impossible to mix up the blasting media

All BEGO 'PLUS' Sandblasters are now equipped with the new push and clean filter system - Filtermodul. BEGO 'Z' Sandblasters for connection to central extraction units are also available. Please contact Metrodent for further details.

10% off rrp\*



EasyBlast



DuoStar



KoroStar



Protempomatic

Take a closer look at the full range of BEGO Sandblasters including the EasyBlast Protempomatic KoroStar and DuoStar at [bego.com](http://bego.com)

EasyBlast £1,129.12\*  
Order ref: 26385  
rrp £1,254.58

DuoStar-Plus £2,476.95\*  
Order ref: 26395  
rrp £2,752.17

KoroStar-Plus £1,932.98\*  
Order ref: 26400  
rrp £2,147.75

Protempomatic-Plus £3,179.12\*  
Order ref: 26390  
rrp £3,532.35

\*Introductory offer ends 31/10/2011

ORDER LINE 01484 466700