



→ VALUE FOR MONEY

→ A TRUE CLASSIC

Metrocyl has been formulated by Metrodent for over 60 years. Over this time it has developed into a denture base acrylic of high quality and consistency.

→ METROCRYL ACRYLICS/HEATCURE

Metrocyl cross-linked acrylics manufactured and tested to BS EN ISO 1567: 2001 Type 1 Class 1. Cadmium Free.

Universal Snap stage 15 minutes, Work time 1 hour at 23°C.
Can be quick cured or overnight.

Rapid Cure Snap stage 30 minutes, Work time 30 minutes at 23°C.
Can be quick cured in 20 minutes.

Metrocyl is available in the following shades. Natural, Standard, Veined, ST Veined, Clear and Translucent.

***When ordering any Metrocyl denture base please state the type and shade required.**

Availability	Unit	Net. (£)
Metrocyl Polymer (Powder)*		
Universal or Rapid - Choose the shade required	1 Kg	£17.95
Universal or Rapid - Choose the shade required	3 Kg	£53.58
Universal or Rapid - Choose the shade required	12 Kg	£146.83
Metrocyl Monomer (Liquid)*		
Universal (for use with pigmented denture base powders)	500 ml	£6.89
Universal (for use with pigmented denture base powders)	1 Litre	£12.65
Universal (for use with pigmented denture base powders)	2.5 Litre	£28.80
Universal (for use with pigmented denture base powders)	2 x 2.5 Litre	£40.22
Metrocyl BULK pack*		
Universal or Rapid Cure - Choose the shade required	12Kg/2 x 2.5 Litre	Normal Price £155.25
Metrocyl bought separately would cost £204.45		Special Offer £140.29

→ METRODENT

believes in giving you value for money and in these difficult economic times we are constantly striving to give you a better deal. We can now offer the Metrocyl Bulk Pack at an even greater saving.

→ NOW ONLY

£140.29

(12 Kg powder/2 x 2.5 litre liquid)*



 **Metrodent** → ORDERLINE 01484 466700



Metrodent Limited Lowergate Works, Lowergate, Paddock, Huddersfield HD3 4EP Metrodent.com

Metrodent believes the prices and information contained in this advert are accurate. However the information and prices are subject to change without notice and should not be construed as a commitment by Metrodent Limited. Metrodent assumes no responsibility for any errors that may appear in this advert. All prices are Net, excluding Vat.