



# THE COMPLETE SOLUTION FOR CEMENTATION

G-CEM LinkForce®  
G-CEM LinkAce®  
& GC FujiCEM® 2





## WELCOME TO THE GC FAMILY OF CEMENTS

There are many options to consider when planning and delivering indirect restorations. The selection of the optimal cementation system is the guarantee for a beautiful smile and long-term patient satisfaction. GC America offers both glass ionomer and composite resin based cement solutions to meet any of your cementation challenges.



### G-CEM LinkForce® Dual-Cure Adhesive Resin Cement

- A strong, esthetic, universal solution for any adhesive cementations
- The ideal cement to use with our **NEW GC Initial® LiSi Press Lithium Disilicate Ingots** and our **CERASMART® Force Absorbing Hybrid CAD/CAM Blocks**
- Four shades available in both cement and try-in paste
- Excellent color stability and easy clean-up



### G-CEM LinkAce® Self-Adhesive Resin Cement

- Innovative initiator system (chemical) offers the optimal self-cure mode
- High bond durability to pretreated Zirconia in ONE step thanks to proprietary phosphate monomers
- Amazing color stability thanks to low water sorption and being HEMA-free
- Small homogeneously distributed particles provide excellent wear resistance
- Easy excess cement removal - only 1-2 second tack cure



### GC FujiCEM® 2 Resin-Modified Glass Ionomer Luting Cement


- Used for **17 million restorations\***
- Excellent handling, consistency and clean-up
- Indicated for myriad of indirect restorations made up of zirconia, metals, PFM and resin
- High fluoride release
- Moisture tolerant
- Virtually no post-operative sensitivity

\*Based on units sold.




# ENGINEERED TO ADDRESS YOUR CEMENTATION CHALLENGES


RESTORATION TYPE	G-CEM LinkForce®	G-CEM LinkAce®	GC FujiCEM® 2
PFM and Metal Crowns, Bridges and Posts			
Metal Alloy, Lithium Disilicate, and Zirconia			
Metal, Ceramic, and Fiber Posts			
Composite			
Feldspathic and Leucite Ceramics			
Hybrid/Nano Ceramics			
PRODUCT CHARACTERISTICS	G-CEM LinkForce®	G-CEM LinkAce®	GC FujiCEM® 2
Multiple Shades			
Fluoride Release			
Hydrophilic/Moisture Tolerant			



Source: GCC R&D.



Source: GCC R&D.



Source: GCC R&D.

## A Complete Cementation System

In a constantly changing world of new and innovative dental materials, it is important now more than ever to choose a cementation system you can trust to securely bond any of them. GC's family of cements gives you the flexibility to provide your patients with a comfortable, confident, and trouble-free smile. For a demonstration or to order, contact your local GC sales representative, call our Customer Care Center at 1-800-323-7063, or visit [www.gcamerica.com](http://www.gcamerica.com).

†RelyX™ Unicem Automix, RelyX™ Ultimate, RelyX™ Luting, RelyX™ Luting Plus and RelyX™ Luting Plus AM are trademarks of 3M. SpeedCEM® and Variolink® II are registered trademarks of Ivoclar Vivadent AG. Maxcem Elite™ is a trademark of Kerr Corporation.



# ORDERING INFORMATION

## G-CEM LinkForce®

### 009541 G-CEM LinkForce® System Kit

Contains: 3 G-CEM LinkForce® cement syringes (1 each of: A2, Translucent and Opaque; 8.7 g/5 mL each), 45 GC Automix Tips, 5 GC Automix Endo Tip mixing tips with intraoral tip, 1 bottle G-Premio BOND™ (5 mL), 1 bottle G-Multi PRIMER™ (5 mL), 1 bottle G-Premio BOND™ Dual Cure Activator (3 mL), 1 syringe GC ETCHANT (3.6 mL; available in System Kit only), 25 Metal dispensing tips + 2 covers (for GC ETCHANT & Try-in Paste), 3 syringes G-CEM LinkForce® Try-in Paste (A2, Translucent and Opaque; 1.5 g each), 50 disposable applicators and 20 disposable dispensing dishes.

### G-CEM LinkForce® Starter Kit

Contains: 1 G-CEM LinkForce® cement syringe (8.7 g/5 mL), 1 bottle G-Premio BOND™ (5 mL), 1 bottle G-Multi PRIMER™ (5 mL), 20 GC Automix Tips, 50 disposable applicators and 20 disposable dispensing dishes.

010118 A2                      010119 Translucent

### G-CEM LinkForce® Cement Refills

Contains: 1 syringe (8.7 g/5 mL) and 20 GC Automix Tips.

010120 A2                      010122 Opaque                      010123 Bleach  
010121 Translucent

### G-CEM LinkForce® Try-In Paste Refills

Contains: 1 syringe (1.5 g/1.2 mL), 3 dispensing tips (metal type) and 1 tip cover.

009548 A2                      009550 Opaque                      009551 Bleach  
009549 Translucent

### 009275 G-Premio BOND™ Kit

Contains: 1 bottle (5 mL), 50 disposable applicators and 20 disposable dispensing dishes.

### 009276 G-Premio BOND™ Unit Dose

Contains: 50 unit doses (0.1 mL) and 50 disposable applicators (fine).

### 009277 G-Premio BOND™ Bottle

Contains: 1 bottle (5 mL).

### 009553 G-Multi PRIMER™ Bottle

Contains: 1 bottle (5 mL).

### 009552 G-Premio BOND™ Dual Cure Activator (DCA)

Contains: 1 bottle (3 mL).

### 002902 GC Automix Tips

Contains: 20 tips.

### 002903 GC Automix Tips for Endo

Contains: 10 tips and 10 extension tips.

### 001926 Flo/LoFlo Metal Dispensing Tips

Contains: 20 Metal Dispensing Tips for G-CEM LinkForce® Try-In Paste or GC Etchant.

## G-CEM LinkAce®

### G-CEM LinkAce® Syringe Packages

Contains: Two 4.6 g (2.7 mL) syringes; 15 regular GC Automix Tips and 5 GC Automix Tips for Endo.

004857 A2                      004859 AO3                      004860 BO1  
004858 Translucent

### 002902 GC Automix Tips

Contains: 20 tips.

### 002903 GC Automix Tips for Endo

Contains: 10 tips and 10 extension tips.

## GC FujiCEM® 2

### 441000 (US) / 451000 (Canada) GC FujiCEM® 2 Automix Package

Contains: Two (7.2 mL; 13.3 g) cartridges and 42 mixing tips.

### 004647 GC FujiCEM® 2 Dispenser

Contains: One dispenser.

### 441001 (US) / 451001 (Canada) GC FujiCEM® 2 Starter Kit

Contains: One (7.2 mL; 13.3 g) cartridge, 20 mixing tips and 1 GC FujiCEM® 2 Dispenser.

### 441002 (US) / 451002 (Canada) GC FujiCEM® 2 Refill

Contains: Two (7.2 mL; 13.3 g) cartridges and a #22 Mixing Pad.

### 431578 GC FujiCEM® 2 Mixing Tips SL

Contains: 20 mixing tips.

### 431577 GC FujiCEM® 2 Mixing Tips Intraoral Tip SL

Contains: 10 mixing tips.

