



New GC FujiVest[®] Platinum II from GC

Better flow • longer working time • increased expansion for a better workflow



Top quality investment
material for the creation of
highly precise castings
with precious,
semi-precious and
Pd base alloys.

GC



A more **exact fit** than ever before

Perfect control and a highly accurate fit – these are the most important requirements for a high-precision casting investment material. With **GC Fujivest Platinum II**, GC Europe has set a new standard for crown and bridge fabrication involving precious, semi-precious and palladium-based alloys.

The impressive features of this carbon-free crown and bridge (C&B) precision investment include excellent working properties and perfect control of expansion – for both quick-heating & step-heating procedures.

Smooth surfaces, putting every fine detail in the spotlight

As a high-precision investment material for crown and bridge work, **GC Fujivest Platinum II** is the material of choice when you are looking for the highest precision. Its perfectly balanced composition contains carefully selected fine-particle fillers for a significant reduction of surface roughness. This enables accurate reproduction of the finest details of the wax pattern and enhanced quality of the final casting.




User-friendliness and peace of mind

GC Fujivest Platinum II is easy to work with. The extended working time of up to 6 minutes and improved fluidity of **GC Fujivest Platinum II** provides extra time for you to ensure that the wax pattern is accurately invested. Devesting has never been easier and faster, leading to a reduced risk of damaged castings.



Great versatility in an all-round material

With an increased expansion range, this top class material has become even more versatile. The investment's stable compressive strength ensures edge stability. All of this adds up to a material that is the perfect partner when you are looking for a solution for castings in complicated cases such as implant suprastructures and long-span bridges. Since this material is manufactured by the inventor of Fujivest® Super, the first ever quick-heating investment, it goes without saying that **GC Fujivest Platinum II** is a superior investment material that can be used in both quick-heating & step-heating procedures.



Reliable and consistent high quality for a perfect fit – every time

Consistent quality from one batch to the next is an absolute must – especially in this day and age, when success revolves around uniform results and at a cost you can afford. Ensuring products comply with specifications is, however, also about careful monitoring. That is why we have watertight quality control – careful checks on each and every batch. Not just once, but several times throughout the overall manufacturing process.

All GC Europe investment materials undergo a variety of inspections, from initial raw materials to the finished product. Plus, each inspection is based on standard benchmarking criteria – before the final batch is released.



GC FujiVest® Platinum II

- Specially developed for high-precision castings of precious, semi-precious and Pd-based alloys
- Extended working time, improved fluidity
- Facilitates a perfect fit for complex cases, such as implant suprastructures and long-span bridges
- Extremely smooth casting surfaces
- Easy to steer and very wide expansion range
- Devests easily, with reduced risk of damaging castings



Packages

GC FujiVest® Platinum II	
900671	6 kg powder (100 x 60 g bags)
900672	6 kg powder (67 x 90 g bags)
900673	6 kg powder (40 x 150 g bags)
900676	900 ml liquid

GC EUROPE N.V.

Head Office
 Researchpark
 Haasrode-Leuven 1240
 Interleuvenlaan 33
 B-3001 Leuven
 Tel. +32.16.74.10.00
 Fax. +32.16.40.48.32
 info@gceurope.com
 http://www.gceurope.com

GC UNITED KINGDOM Ltd.

12-15, Coopers Court
 Newport Pagnell
 UK-Bucks. MK16 8JS
 Tel. +44.1908.218.999
 Fax. +44.1908.218.900
 info@uk.gceurope.com
 http://uk.gceurope.com

