

New



Product Guide

Metrodent Exports +44 (0) 1484 466706





Made by Metrodent

Metrocryn Denture Base Acrylics


Over decades Metrodent have developed a range of denture base acrylics for both our UK and Export markets. Metrocryn Rapid Cure Denture Base is available in 9 shade options including clear and ethnic shades. In addition we have an orthodontic resin, Orthomet, for the “spray” technique available in clear and tinted pink.

All our heat cure resins conform to Standard BS EN ISO20795-1: 2008

Metrocryn - Rapid Cure

Denture Base Acrylics

ST Veined (Pink Translucent Veined)



Our most popular acrylic formulated in response to market demand to match the leading brand product. A silky dough is produced with exceptionally good flow properties resulting in reduced flash thickness.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 230C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Stock Code	Description
Ref: R25SV	Metrocryn - Rapid Cure ST Veined 25 Kg
Ref: R50SV	Metrocryn - Rapid Cure ST Veined 50 Kg
Ref: LC25	Metrocryn Clear Liquid 25 Litres
Ref: LC180	Metrocryn Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

ST Veined

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**




Stock Code	Description
Ref: MP25SV	Metropair ST Veined Powder 25 Kg
Ref: MP50SV	Metropair ST Veined Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocryn - Rapid Cure

Denture Base Acrylics

Translucent Veined (Medium Pink Veined)



A popular translucent denture base with veining. Slightly whiter or more opaque than either our Translucent or Light Veined acrylic.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23oC
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Stock Code	Description
Ref: R25TV	Metrocryn - Rapid Cure TV Translucent Veined 25 Kg
Ref: R50TV	Metrocryn - Rapid Cure TV Translucent Veined 50 Kg
Ref: LC25	Metrocryn Clear Liquid 25 Litres
Ref: LC180	Metrocryn Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

TV Translucent Veined

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**




Stock Code	Description
Ref: MP25TV	Metropair TV Translucent Veined Powder 25 Kg
Ref: MP50TV	Metropair TV Translucent Veined Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocyl - Rapid Cure

Denture Base Acrylics

Standard Pink (Natural Pink)



With greater opacity this acrylic has no veining. A lustrous natural pink mimicking healthy gingival mucosa, un-veined.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Stock Code	Description
Ref: R25S	Metrocyl - Rapid Cure Standard Pink 25 Kg
Ref: R50S	Metrocyl - Rapid Cure Standard Pink 50 Kg
Ref: LC25	Metrocyl Clear Liquid 25 Litres
Ref: LC180	Metrocyl Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

Standard Pink

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**




Stock Code	Description
Ref: MP25S	Metropair Standard Pink Powder 25 Kg
Ref: MP50S	Metropair Standard Pink Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocyl - Rapid Cure

Denture Base Acrylics

Translucent Non Veined (Light Gingival Pink)



Translucent denture base, translucency allows the patient's own gingival tones to "colour" the denture base.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Stock Code	Description
Ref: R25T	Metrocyl - Rapid Cure Translucent Non Veined 25 Kg
Ref: R50T	Metrocyl - Rapid Cure Translucent Non Veined 50 Kg
Ref: LC25	Metrocyl Clear Liquid 25 Litres
Ref: LC180	Metrocyl Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

Translucent Non Veined

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**




Stock Code	Description
Ref: MP25T	Metropair Translucent Non Veined Powder 25 Kg
Ref: MP50T	Metropair Translucent Non Veined Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocyl - Rapid Cure

Denture Base Acrylics

Natural (Medium Pink)



A pink opaque denture base mimicking certain gingival tissues. Medium pink with a very subtle orange tint and good opacity.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Stock Code	Description
Ref: R25N	Metrocyl - Rapid Cure Natural 25 Kg
Ref: R50N	Metrocyl - Rapid Cure Natural 50 Kg
Ref: LC25	Metrocyl Clear Liquid 25 Litres
Ref: LC180	Metrocyl Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

Natural

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**



Stock Code	Description
Ref: MP25N	Metropair Natural Powder 25 Kg
Ref: MP50N	Metropair Natural Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocyl - Rapid Cure

Denture Base Acrylics

Ethnic (Dark Purple shade of Pink veined)



Developed to match the gingival tissues of patients with afro-caribbean heritage.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Made by Metrodent



Stock Code	Description
Ref: R25E	Metrocyl - Rapid Cure Ethnic 25 Kg
Ref: R50E	Metrocyl - Rapid Cure Ethnic 50 Kg
Ref: LC25	Metrocyl Clear Liquid 25 Litres
Ref: LC180	Metrocyl Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

Ethnic

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**



Made by Metrodent




Stock Code	Description
Ref: MP25E	Metropair Ethnic Powder 25 Kg
Ref: MP50E	Metropair Ethnic Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocyl - Rapid Cure

Denture Base Acrylics

Dark V1 Veined (A more reddish Pink Veined)



Developed originally for export to our East Asian markets and offering a variation in gingival tissue matching.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Stock Code	Description
Ref: R25DV	Metrocyl - Rapid Cure Dark V1 Veined 25 Kg
Ref: R50DV	Metrocyl - Rapid Cure Dark V1 Veined 50 Kg
Ref: LC25	Metrocyl Clear Liquid 25 Litres
Ref: LC180	Metrocyl Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

Dark V1 Veined

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**




Stock Code	Description
Ref: MP25DV	Metropair Dark V1 Veined Powder 25 Kg
Ref: MP50DV	Metropair Dark V1 Veined Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocyl - Rapid Cure

Denture Base Acrylics

Light Veined (Light Translucent Pink Veined)



Similar to our most popular ST Veined acrylic but in a lighter more subtle pink suitable for patients with paler gingival tissues.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Made by Metrodent



Stock Code	Description
Ref: R25LV	Metrocyl - Rapid Cure Light Veined 25 Kg
Ref: R50LV	Metrocyl - Rapid Cure Light Veined 50 Kg
Ref: LC25	Metrocyl Clear Liquid 25 Litres
Ref: LC180	Metrocyl Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

Light Veined

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**



Made by Metrodent



Stock Code	Description
Ref: MP25LV	Metropair Light Veined Powder 25 Kg
Ref: MP50LV	Metropair Light Veined Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Metrocryn - Rapid Cure

Denture Base Acrylics

Clear

Ideal for clear palates, TMJ appliances, orthodontic appliances, surgical guides etc. Using rapid cure gives an exceptionally clear acrylic whereas universal can sometimes have a slight opacity.

Advantages

- ✓ Snap stage 30 minutes
- ✓ Work time 30 minutes at 23°C
- ✓ Can be quick cured in 20 minutes

Universal Liquid

- ✓ For all coloured acrylics



Made by Metrodent



Stock Code	Description
Ref: R25C	Metrocryn - Rapid Cure Clear 25 Kg
Ref: R50C	Metrocryn - Rapid Cure Clear 50 Kg
Ref: LC25	Metrocryn Clear Liquid 25 Litres
Ref: LC180	Metrocryn Clear Liquid 180 Litres

Metropair

Denture Repair Acrylic

Clear

Cold cure acrylic for denture repairs, additions, relines etc., non slumping, quick setting and free of porosity if cured in a pressure flask. **Use only with Metropair Liquid.**



Made by Metrodent



Stock Code	Description
Ref: MP25C	Metropair Clear Powder 25 Kg
Ref: MP50C	Metropair Clear Powder 50 Kg
Ref: MPL25	Metropair Liquid 25 Litres
Ref: MPL180	Metropair Liquid 180 Litres

Orthomet

Orthodontic Acrylic Clear

Orthomet is used in the “spray” technique for the manufacture of Orthodontic appliances. With the use of a pressure flask such as the Metroclave or Hi-Lo the material will be cured in approximately 10 minutes showing no porosity or shrinkage.



Stock Code	Description
Ref: OP25	Orthomet Powder 25 Kg
Ref: OP50	Orthomet Powder 50 Kg
Ref: OLC25	Orthomet Clear Liquid 25 Litres



Custom Impression Tray Acrylic

A high quality self cure material with excellent working properties and economic in use.

Blue - standard time.
Dough time 2 to 3 minutes.

Green - long work time.
Dough time 4 to 5 minutes..



Stock Code	Description
Ref: MTFB25	Metroform Blue Powder 25 Kg
Ref: MTFB50	Metroform Blue Powder 50 Kg
Ref: MTFG25	Metroform Green Powder 25 Kg
Ref: MTFG50	Metroform Green Powder 50 Kg
Ref: MTFL25	Metroform Liquid (for the above powders) 25 Kg

PCS (Plaster Coating Solution)



Cold Mould Seal

Premium grade formaldehyde free alginate solution for insulation of denture base resins from gypsum moulds and casts when making dentures, acrylic crowns etc.



Stock Code	Description
Ref: PCS1	PCS (Plaster Coating Solution) 1 Litre
Ref: PCS2	PCS (Plaster Coating Solution) 2 Litres
Ref: PCS5	PCS (Plaster Coating Solution) 5 Litres
Ref: PCS25	PCS (Plaster Coating Solution) 25 Litres
Ref: PCS200	PCS (Plaster Coating Solution) 200 Litres



Duplicating Gels

Reversible Hydrocolloid



Plastergel

Plaster compatible and gives a very smooth investment of plaster model surface.
Available in Red or Clear



Metrogel

Proven reversible hydrocolloid gel with good elasticity. Ideal for partial denture work.
Available in Blue or Clear



Supergel

Tougher than Metrogel means that it is particularly suitable for attachment work and demanding partial denture work.
Available in Green or Clear



Stock Code	Description
Ref: PGEL	Plastergel - Duplicating Gel, Red - 5 Kg
Ref: PGEL/C	Plastergel - Duplicating Gel, Clear - 5 Kg
Ref: MGEL	Metrogel - Duplicating Gel, Blue - 5 Kg
Ref: MGEL/C	Metrogel - Duplicating Gel, Clear - 5 Kg
Ref: SGEL	Supergel - Duplicating Gel, Green - 5 Kg
Ref: SGEL/C	Supergel - Duplicating Gel, Clear - 5 Kg

For the manufacture of acrylic teeth. Kallodoc is available in basic colours and pre-blended shades



Stock Code	Description
Ref: KAL50/000	Kallodoc Acrylic Tooth Powder Basic Colour - Clear 000 - 50 Kg
Ref: KAL50/001	Kallodoc Acrylic Tooth Powder Basic Colour - Opal 001 - 50 Kg
Ref: KAL50/103	Kallodoc Acrylic Tooth Powder Basic Colour - Cream 103 - 50 Kg
Ref: KAL50/202	Kallodoc Acrylic Tooth Powder Basic Colour - Yellow 202 - 50 Kg
Ref: KAL50/453	Kallodoc Acrylic Tooth Powder Basic Colour - Pink 453 - 50 Kg
Ref: KAL50/502	Kallodoc Acrylic Tooth Powder Basic Colour - Brown 502 - 50 Kg
Ref: KAL50/703	Kallodoc Acrylic Tooth Powder Basic Colour - Blue 703 - 50 Kg
Ref: KAL50/952	Kallodoc Acrylic Tooth Powder Basic Colour - Grey 952 - 50 Kg
Ref: KAL50/170	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Incisal 170 - 50 Kg
Ref: KAL50/171	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Cream 171 - 50 Kg
Ref: KAL50/173	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Light Ivory 173 - 50 Kg
Ref: KAL50/174	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Light Ivory 174 - 50 Kg
Ref: KAL50/175	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Medium Ivory 175 - 50 Kg
Ref: KAL50/177	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Medium Ivory 177 - 50 Kg
Ref: KAL50/178	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Dark Ivory 178 - 50 Kg
Ref: KAL50/179	Kallodoc Acrylic Tooth Powder Pre Blended Shade - Dark Ivory 179 - 50 Kg




Made by Metrodent

Modelling Waxes


Metrodent have been a manufacturer of modelling waxes since 1936 and produce a comprehensive range used and relied upon by Dental Technicians both in the UK and throughout the world. Using the finest raw materials each of our waxes is blended and produced to strict quality control standards giving a first class British made product.

Modelling Waxes




Metrowax No.1

Our highest melting point modelling wax




Metrowax No.2

Popular blend




Metrowax No.4

No.2 formula but richer in colour




Metrowax No.6

Blended with natural and synthetic waxes




Metrowax No.8

Blended with natural and synthetic waxes



Light Pink

No.2 formula but lighter




Metromorphic

Additional Bees Wax creates a pleasant aroma



Stock Code	Description
Ref: 1WAX500	Metrowax No.1 - 500g
Ref: 1WAX2.5	Metrowax No.1 - 2.5 Kg
Ref: 1WAX20	Metrowax No.1 - 20 Kg
Ref: 2WAX500	Metrowax No.2 - 500g
Ref: 2WAX2.5	Metrowax No.2 - 2.5 Kg
Ref: 2WAX20	Metrowax No.2 - 20 Kg
Ref: 4WAX500	Metrowax No.4 - 500g
Ref: 4WAX2.5	Metrowax No.4 - 2.5 Kg
Ref: 4WAX20	Metrowax No.4 - 20 Kg
Ref: 6WAX500	Metrowax No.6 - 500g
Ref: 6WAX2.5	Metrowax No.6 - 2.5 Kg
Ref: 6WAX20	Metrowax No.6 - 20 Kg
Ref: 8WAX500	Metrowax No.8 - 500g
Ref: 8WAX2.5	Metrowax No.8 - 2.5 Kg
Ref: 8WAX20	Metrowax No.8 - 20 Kg
Ref: LPWAX500	Light Pink - 500g
Ref: LPWAX2.5	Light Pink - 2.5 Kg
Ref: LPWAX20	Light Pink - 20 Kg
Ref: 7WAX500	Metromorphic - 500g
Ref: 7WAX2.5	Metromorphic - 2.5 Kg
Ref: 7WAX20	Metromorphic - 20 Kg

Modelling Waxes




College Wax

Produced with economy in mind, and in contrast to our virgin melt range, college wax fulfils the needs of a market where price is of prime concern. Ideal for Dental Schools, Colleges and Universities this wax is made of secondary blends of our primary waxes and therefore may have slight variations in colour and handling properties. College wax represents excellent value for money.

Sheet dimensions (Approx) 183 x 87 x 1.75 mm



Stock Code	Description
Ref: CWAX	College Wax - 500g
Ref: CWAX20	College Wax - 20 Kg



Gnathic - Modelling Waxes

Ivory wax sheets and preformed ivory wax bite blocks for superior occlusal analysis.

- > Differentiate your private work at the bite rim stage
- > More accurately represents the dentate area of the rim
- > Improved patient communication
- > Intra oral contrast makes the registration process easier



Stock Code	Description
Ref: GWAX500	College Wax - 500g
Ref: GWAXBB48	College Wax - 20 Kg

Wax Bite Blocks

For detailed occlusal records

Block dimensions 13 x 10 x 114 mm
48 blocks per box



Stock Code	Description
Ref: WAXBB48	Wax Bite Blocks

Ribbon Wax

For dorsal damming and other intra oral applications

126 Strips per box (400 g)
Dimensions 8 x 1.5 x 210 mm



Stock Code	Description
Ref: RWAX	Ribbon Wax - 400g

Wax-Up Sticks

For waxing and setting up try-in dentures

64 Strips per box (500 g)
Dimensions ø 9 x 110 mm



Stock Code	Description
Ref: WUS	Wax-Up Sticks - 500g

Strip Wax

Alternative to Bite Registration Wax

35 Strips per box (200 g)
Dimensions ø 9 x 110 mm



Stock Code	Description
Ref: SWAX	Strip Wax - 200g

Soft Bite Wax

Economical Bite Registration Strips

24 Sticks per box.
Dimensions 12 x 7.5 x 24 mm



Stock Code	Description
Ref: SRBW	Soft Bite Wax

Boxing-In Wax

For boxing impressions prior to casting

18 per box.
Dimensions 300 x 40 x 2.5 mm



Stock Code	Description
Ref: BIW	Boxing-In Wax

Sticky (Cement) Wax



For plaster, metals and resins

Ideal for sticking to plaster, metals and resin.
Dimensions ø 9 x 102 mm



Stock Code	Description
Ref: SW70G	Sticky (Cement) Wax 70g
Ref: SW500G	Sticky (Cement) Wax 500g

Carving Wax Blocks



Used as a training aid by students

For the study of morphology



Stock Code	Description
Ref: DWS511	Carving Wax Block Wax 30 pieces

Dipping Wax



Granular

Cost effective way of producing crown and bridge copings.



Stock Code	Description
Ref: MDWG1KG	Dipping Wax - 1 Kg

Waxes for Modelling



Cervical Wax - Red

45g tin

Crown Wax - Blue

45g tin

Undercut Wax - Pink

45g tin



Stock Code	Description
Ref: W24350	Cervical Wax - 45g
Ref: W24511	Crown Wax - 45g
Ref: W24540	Undercut Wax - 45g

Wax Wire for Sprues



Easy to shape and burn-out without residue

Green 2/3 hard, 250g spool



Stock Code	Description
Ref: W12120	Wax Wire for Sprues - ø 2 mm
Ref: W12125	Wax Wire for Sprues - ø 2.5 mm
Ref: W12130	Wax Wire for Sprues - ø 3 mm
Ref: W12135	Wax Wire for Sprues - ø 3.5 mm
Ref: W12140	Wax Wire for Sprues - ø 4 mm
Ref: W12150	Wax Wire for Sprues - ø 5 mm



Wax Grid Retention Diagonal



For maxillary partial denture frames

Wax grid retentions permit the simple and effective shaping of retentions to total or partial dentures. They guarantee a high level of security in the connection between the resin and the partial denture plate. The large grid retentions facilitate very economical use of material.

Colour: Red, 10 sheets
Dimensions 150 x 75 mm



Stock Code	Description
Ref: W13030	Wax Grid Retention - Diagonal

Wax Grid Retention with holes



For upper jaw

Can be used as retentions for partial maxillary dentures and as a reinforcement for acrylic full maxillary acrylic dentures.

Colour: Red, 25 sheets
Dimensions 65 x 75 mm



Stock Code	Description
Ref: W13020	Wax Grid Retention with holes

Wax Grid Retention Small Mesh



For maxillary partial denture frames

Wax grid retentions permit the simple and effective shaping of retentions to total or partial dentures. They guarantee a high level of security in the connection between the resin and the partial denture plate. The large grid retentions facilitate very economical use of material.

Colour: Red, 25 sheets
Dimensions 40 x 60 mm



Stock Code	Description
Ref: W13000	Wax Grid Retention - Small Mesh

Wax Grid Retention Large Mesh



For maxillary partial denture frames

Wax grid retentions permit the simple and effective shaping of retentions to total or partial dentures. They guarantee a high level of security in the connection between the resin and the partial denture plate. The large grid retentions facilitate very economical use of material.

Colour: Red, 25 sheets
Dimensions 100 x 100 mm



Stock Code	Description
Ref: W13010	Wax Grid Retention - Large Mesh

Wax Border Strips Small with retentions



Time saving

Time savings when modelling upper-jaw frames with large saddles. A great advantage is that the border strip can easily be shaped as desired since the size can be varied by trimming the tips of the retentions. The wax is so supple that it can be shaped easily and reliably as required onto the duplicate model.

Colour: Red, 15 pieces (100 mm length)



Stock Code	Description
Ref: W14010	Wax Border Strip - Small with retentions

Wax Border Strips Large with retentions



Time saving

Time savings when modelling upper-jaw frames with large saddles. A great advantage is that the border strip can easily be shaped as desired since the size can be varied by trimming the tips of the retentions. The wax is so supple that it can be shaped easily and reliably as required onto the duplicate model.

Colour: Red, 15 pieces (130 mm length)



Stock Code	Description
Ref: W14020	Wax Border Strip - Large with retentions

Wax Edge Strips with retentions



Time saving

For the secure attachment of plastic saddles to lower partial dentures.

Colour: Red, 25 pieces (120 mm length)



Stock Code	Description
Ref: W14110	Wax Edge Strips with retentions

Wax Edge Strips with perforated retentions



For lower jaw partial denture frames

For the secure attachment of plastic saddles to lower partial dentures.

Colour: Red, 25 pieces (125 mm length)



Stock Code	Description
Ref: W14120	Wax Edge Strips with perforated retentions

Wax Edge Strips with round holes



For lower jaw partial denture frames

For the secure attachment of plastic saddles to lower partial dentures.

Colour: Red, 25 pieces (125 mm length)



Stock Code	Description
Ref: W13040	Wax Edge Strips with round holes

Wax Profile Sticks Half Round (Continuous Clasps)



Extremely versatile

Extremely versatile multiple applications. Can be used for secondary spruing to clasps, forming continuous clasps or minor connectors.

Colour: Green 50g
Dimensions 1.8 x 0.9 x 145 mm
Dimensions 2.0 x 1.0 x 145 mm



Stock Code	Description
Ref: W12212	Wax Profile Sticks - Half Round (continuous clasps) 1.8 mm
Ref: W12213	Wax Profile Sticks - Half Round (continuous clasps) 2.0 mm

Preparation Wax



For the partial denture technique

The preparation wax is exceptionally malleable, allowing it to be adapted to the master model perfectly and with firm adhesion, which saves having to use an additional wax adhesive. The exemplary shape retention and edge strength of the preparation wax, with a high solidification point of approx. 70°C, mean that it can be used with duplicating hydrocolloid at working temperatures of 55°C. Simple removal from the master model following duplication rounds off the userfriendly working characteristics perfectly.

Colour: Red (15 sheets)
Dimensions 175 x 80 mm



Stock Code	Description
Ref: W11750	Preparation Wax 0.5 mm
Ref: W11760	Preparation Wax 0.6 mm
Ref: W11770	Preparation Wax 0.7 mm

Smooth Casting Wax



Simple crease free adaptation

Simple, crease-free adaptation. Adheres firmly to the investment model and burns out leaving no residue. The high transparency of the wax makes for optimal clarity of the construction markings on the master model and saves unnecessary, time-consuming corrections to the wax-up.

Colour: Green (15 sheets)
Dimensions 175 x 80 mm



Stock Code	Description
Ref: W10030	Smooth Casting Wax 0.30 mm
Ref: W10040	Smooth Casting Wax 0.40 mm
Ref: W10050	Smooth Casting Wax 0.50 mm
Ref: W10060	Smooth Casting Wax 0.60 mm
Ref: W10070	Smooth Casting Wax 0.70 mm

Stippled Casting Wax



Tried and tested

Tried and tested wax for modelling the bases of upper partial dentures. Can be easily adapted and adheres firmly to the investment model with no additional wax adhesive
The stippled casting wax is available in three different surface textures – from fine to coarse – and allows customisation of the surface shape as required by the practitioner.
The individual stippling of the cast partial denture base facilitates the gripping of food and reduces the foreign body sensation for the patient's tongue.



Colour: Green (15 sheets)
Dimensions 150 x 75 mm

Stock Code	Description
Ref: W11035	Stippled Casting Wax - Fine 0.35 mm
Ref: W11040	Stippled Casting Wax - Fine 0.40 mm
Ref: W11050	Stippled Casting Wax - Fine 0.50 mm
Ref: W11060	Stippled Casting Wax - Fine 0.60 mm
Ref: W11340	Stippled Casting Wax - Medium 0.40 mm
Ref: W11140	Stippled Casting Wax - Coarse 0.40 mm
Ref: W11150	Stippled Casting Wax - Coarse 0.50 mm
Ref: W11160	Stippled Casting Wax - Coarse 0.60 mm



Wax Profile Sticks Lower Jaw Bars



The anatomical design of these profiles ensures excellent patient acceptance with the finished prosthesis

Colour: Green, (65g)
Dimensions 145 mm length



Stock Code	Description
Ref: W12250	Wax Profile Sticks Lower Jaw Bars 4.0 x 2.0 mm
Ref: W12251	Wax Profile Sticks Lower Jaw Bars 4.0 x 1.75 mm
Ref: W12252	Wax Profile Sticks Lower Jaw Bars 4.3 x 2.3 mm
Ref: W12240	Wax Profile Sticks Upper Jaw Bars 5.0 x 1.3 mm

Wax Patterns/Wax Clasp Profiles



For Bonyhard Type A

These preshaped clasp profiles simplify modelling and help save time. The shape of the profiles enables a large number of variations by shortening the wax form.

Colour: Green
(10 sheets) 200 clasps



Stock Code	Description
Ref: W50020	Wax Patterns/Wax Clasp Profiles for Bonyhard Type A

Wax Patterns/Wax Clasp Profiles



For Bonyhard Type B

These preshaped clasp profiles simplify modelling and help save time. The shape of the profiles enables a large number of variations by shortening the wax form.

Colour: Green
(10 sheets) 200 clasps



Stock Code	Description
Ref: W50021	Wax Patterns/Wax Clasp Profiles for Bonyhard Type B

Wax Patterns/Wax Clasp Profiles



For Premolars

These preshaped clasp profiles simplify modelling and help save time. The shape of the profiles enables a large number of variations by shortening the wax form.

Colour: Green
(10 sheets) 200 clasps



Stock Code	Description
Ref: W50010	Wax Patterns/Wax Clasp Profiles for Premolars

Wax Patterns/Wax Clasp Profiles



For Molars

These preshaped clasp profiles simplify modelling and help save time. The shape of the profiles enables a large number of variations by shortening the wax form.

Colour: Green
(10 sheets) 200 clasps



Stock Code	Description
Ref: W50000	Wax Patterns/Wax Clasp Profiles for Molars

Wax Patterns/Ring Clasp Profile - Straight



Slender standard clasp profiles for molars.

Colour: Green
(10 sheets) 200 clasps



Stock Code	Description
Ref: W50031	Wax Patterns/Ring Clasp Profile - Straight

Wax Patterns/Ring Clasp Profile - Curved



Slender standard clasp profiles for molars.

Colour: Green
(10 sheets) 200 clasps



Stock Code	Description
Ref: W50030	Wax Patterns/Ring Clasp Profile - Curved

Wax Profile Sticks - Beading Wires



The malleable nature of these profiles ensure a neat and accurate finish to saddle areas.

Colour: Green, 145 mm length



Stock Code	Description
Ref: W12201	Wax Profile Sticks - Beading Wire 25g ø 0.8 mm
Ref: W12203	Wax Profile Sticks - Beading Wire 25g ø 1.0 mm
Ref: W12207	Wax Profile Sticks - Beading Wire 25g ø 1.5 mm
Ref: W12208	Wax Profile Sticks - Beading Wire 25g ø 2.0 mm

Metrodent Casting Waxes

Partial Upper



Wax Patterns/ Wax Clasp Profile for Bonyhard Type A
Ref: W50020



Wax Profile Sticks - Half Round (Continuous Clasps)
Ref: W12212




Wax Profile Sticks Beading Wire - 1 mm
Ref: W12203



Wax Grid Retention with holes for upper jaw
Ref: W13020



Wax Patterns/Ring Clasp Profile - Straight
Ref: W50031



Wax Patterns/ Wax Clasp Profile for molars
Ref: W50000




Stippled Casting Wax - Fine 0.4 mm
Ref: W11040



Flat Casting Strips Upper Jaw
Ref: 40461

Metrodent Casting Waxes

Lower Single Bar




Metrodent Casting Waxes


Full Upper



Wax Border Strips - Large
with retentions
Ref: W14020



Wax Edge Strips with retentions
Ref: W14110



Stippled Casting Wax - Fine 0.4 mm
Ref: W11040



Flat Casting Strips
Upper Jaw
Ref: 40461

Denture Boxes

For the storage of mouthguards, dentures and retainers. Approximately 45 mm deep.



Grey



Light Blue



Blue



Pink



Vanilla



White



Rose



Yellow

Stock Code	Description
Ref: M13M-XGY	Grey Denture Box - Pack 10
Ref: M13M-XOB	Light Blue Denture Box - Pack 10
Ref: M13M-XDB	Blue Denture Box - Pack 10
Ref: M13M-XPk	Pink Denture Box - Pack 10
Ref: M13M-XFV	Vanilla Denture Box - Pack 10
Ref: M13M-XWH	White Denture Box - Pack 10
Ref: M13M-XDR	Rose Denture Box - Pack 10
Ref: M13M-XYE	Yellow Denture Box - Pack 10

Ortho Boxes

For the storage of orthodontic appliances.

Approximately 20 mm deep.



Dark Blue



Marble Pink



Light Blue



White



Green



Marble Grey



Red



Yellow




Beige

Stock Code	Description
Ref: M131-XDB	Dark Blue Ortho Box - Pack 10
Ref: M131-MBK/WH	Marble Pink Ortho Box - Pack 10
Ref: M131-XBL	Light Blue Ortho Box - Pack 10
Ref: M131-XWH	White Ortho Box - Pack 10
Ref: M131-XGR	Green Ortho Box - Pack 10
Ref: M131-XPR	Marble Grey Ortho Box - Pack 10
Ref: M131-XRD	Red Ortho Box - Pack 10
Ref: M131-XYE	Yellow Ortho Box - Pack 10
Ref: M131-XBE	Beige Ortho Box - Pack 10


Dento/Ortho Boxes

For the storage of prosthetic and orthodontic appliances.
High quality robust plastic hinged containers.

Approximately 40 mm deep.



Blue



Dark Blue



White



Red




Pink



Green



Bright Green



Yellow

Stock Code	Description
Ref: HW355142	Blue Dento/Ortho Box - Pack 10
Ref: GLG3BM	Dark Blue Dento/Ortho Box - Pack 10
Ref: GLG3B	White Dento/Ortho Box - Pack 10
Ref: GLG3RF	Red Dento/Ortho Box - Pack 10
Ref: HW355143	Pink Dento/Ortho Box - Pack 10
Ref: GLG3V	Green Dento/Ortho Box - Pack 10
Ref: GLG3VM	Bright Green Dento/Ortho Box - Pack 10
Ref: GLG3G	Yellow Dento/Ortho Box - Pack 10

Despatch Boxes

Available in small and large sizes.
Each box comes with 2 pre cut foam inserts.

Small Dimensions 29 x 29 x 22 mm
Large Dimensions 55 x 45 x 22 mm



Stock Code	Description
Ref: KEY1160010	Small Despatch Box - Pack 10
Ref: KEY1160020	Large Despatch Box - Pack 10

Polishing Trough

Ideal for wet pumice polishing (without suction port).
Constructed of hardwearig rubber.



Stock Code	Description
Ref: DPT	Polishing Trough (Rubber)

Casting Rings and Casting Ring Bases

Rubber bases and stainless steel rings available in 3 sizes.



Stock Code	Description
Ref: P7110	Casting Ring Base - No.1
Ref: P7140	Casting Ring Base - No.3
Ref: P7150	Casting Ring Base - No.6
Ref: SSCR1	Casting Ring ø 30 x 40 mm - No.1
Ref: SSCR3	Casting Ring ø 48 x 55 mm - No.3
Ref: SSCR6	Casting Ring ø 65 x 55 mm - No.6

Glass Mortar and Pestle

Grounded and etched glass.
2 inch diameter



Stock Code	Description
Ref: P5035	Glass Mortar and Pestle Set

Cotton Wool Holder


Cotton wool holder with spring loading platform



Stock Code	Description
Ref: P5036	Cotton Wool Holder

Maun Cutters

Compound action, hared steel nose



Stock Code	Description
Ref: MC	Maun Cutters



Dowel Pins

No.2 medium standard flat sided with knurled base
1000 Dowel pins per box



Stock Code	Description
Ref: DPIN2	No.2 Dowel Pins x 1000

Cold Cure Vessels

The simple and reliable solution

Metrodent have been making the Metroclave and Hi-Lo Cold Cure Vessels for over 30 years. We think they are the best, and many Technicians agree with us. Part of their success is based on their simplicity. To process cold cure resins two things are required heat and pressure. Using the Metroclave and Hi-Lo the heat is supplied by filling the vessels with warm water and th pressure is applied by screwing down the lid and indicated on the gauge. That's it!

Advantages

- ✓ No compressor required
- ✓ No heating plate needed
- ✓ Does not tie-up the flask press whilst in use
- ✓ Requires little maintenance
- ✓ Gives years of service
- ✓ Inexpensive

Metroclave

The lid is recessed in a well in the body of the vessel, this allows it to be used on a bench away from the sink or drainer. If the Metroclave is filled to just the right level, the well allows the lid to be fitted without spillage. Curing time 6 to 15 minutes according to the water temperature. Holds up to 5 models. Dimensions: ø 160 x 90 mm.

Hi-Lo

This is smaller of the two vessels and has a capacity for up to three models. It is ideal for the smaller laboratory or department or even for occasional use in a dental surgery. Curing time 6 to 15 minutes according to the water temperature. Dimensions; ø 110 x 110 mm.



Stock Code	Description
Ref: MTC	Metroclave Cold Cure Vessel
Ref: HLA	Hi-Lo Cold Cure Vessel

Universal Type A Articulator

Ideal for easy positioning of models without the split cast system.



Stock Code	Description
Ref: P5150	Universal Type A Articulator

Free Plane Articulator

Lightweight anatomic

Made of diecast aluminium. It can be used following the geometrical laws of medium values (Bonwill triangle). Fixing according to centric. Automatic return to the hinge axis with springs. Perfect stability also in reversed position. Easy access and visibility to the lingual surfaces. Interchangeable self centring aluminium mounting plates (ref 5000-1) or magnetic plates (ref 5000-1M). Pure translation movement (key and instructions included). Anti-shock hinge rod in case articulator falls accidentally. Bigger size lock screws for mounting plates.



Stock Code	Description
Ref: ATA	Free Plane Articulator
Ref: S5000-2	Occlusal table
Ref: S5000-52	Guide rod lock screws
Ref: ATAPN	Mounting plate screw (S-50006)
Ref: S5000-8	Incisal indicator pin

Galetti® Articulator

Plasterless fixed plane articulator


The worlds most recognised Plaster - Free Articulator. Solidly constructed from Aluminium and Brass. Hard wearing and corrosion resistant. Any type of dental model can be placed onto the Galetti articulator without the application of mounting plaster.

This enables precise and controlled adjustment of the occlusion to exacting tolerance. Both upper and lower jaws may be separated allowing for greater access and stability.

The Galetti articulator enables the user to check occlusal and articulated balance for several types of dental prostheses. Constant assessment of mandibular motion may be verified by adjustment of the articulator.

Advantages

- ✓ Check functional bite registration
- ✓ Assess Inter Cuspal Position (Centric Occlusion)
- ✓ Assist with planning occlusal schemes
- ✓ Adjust occlusal interferences of dental prostheses



Stock Code	Description
Ref: ATE	Galetti plasterless fixed plane articulator
Ref: K60492	1. Later claw bolt
Ref: K60484	2. Moveable locking claw bolt
Ref: K60481	3. Tension bolt
Ref: K60487	4. Adjusting bolt
Ref: K60488	5. Knurled lock nut
Ref: K60489	6. Knurled nut
Ref: K60496	7. Tension spindle
Ref: K60494	8. Bolt for wing nut
Ref: K60479	9. Adjusting screw
Ref: K60480	10. Locking device
Ref: K60495	11. Ball shaft

Articulating Paper

12 Booklets
Colour: Blue



Stock Code	Description
Ref: SUN102	Articulating Paper (Blue)

Z.O.E. Impression Paste

Kit

Contents:
Base - Zinc Oxide Paste - White 150g
Activator - Eugenol Paste - Red 60g



Stock Code	Description
Ref: SSWIP	Z.O.E. Impression Paste Kit

Denture ID Kit



1 metre x radiopaque metal strip and permanent marker pen.
Includes instructions. Embedded into the palatal or lingual
aspect of a denture for individual Identification.



Stock Code	Description
Ref: DIK	Denture ID Kit

Tracing Sticks



Impression Compound

For adjusting the periphery of impression trays.

Colour: Green
Dimensions ø8 x 115 mm



Stock Code	Description
Ref: COM621	Tracing Sticks (15 per box)

Wapo Bars

Strengthening bars for incorporating in acrylic dentures



Stock Code	Description
Ref: WBL/6	Wapo bar lower - 6 cm
Ref: WBL/7	Wapo bar lower - 7 cm
Ref: WBL/8	Wapo bar lower - 8 cm
Ref: WBU/7	Wapo bar upper - 7 cm

Lingual Bars (Preformed)

Preformed lingual bars available in three different lengths, offering the dental technician significant time saving.

Pack 10



Stock Code	Description
Ref: LB6	Small - 6 cm (10 Pack)
Ref: LB7	Medium - 7 cm (10 Pack)
Ref: LB8	Large - 8 cm (10 Pack)

Ball End Clasps

Alternative forms of clasp retention for partial dentures and orthodontic appliances.

Pack 20



Stock Code	Description
Ref: BEC7	Ball End Clasps - 0.7 mm (20 Pack)
Ref: BEC8	Ball End Clasps - 0.8 mm (20 Pack)
Ref: BEC9	Ball End Clasps - 0.9 mm (20 Pack)

Roach T and J Clasps

Alternative forms of clasp retention for partial dentures and orthodontic appliances.

Roach T Clasps Pack 10
J-Shaped Clasps Pack 20



Stock Code	Description
Ref: RCS	Roach T Clasp - Small (20 Pack)
Ref: RCM	Roach T Clasp - Medium (20 Pack)
Ref: RCL	Roach T Clasp - Large (20 Pack)
Ref: JC8	J Clasp - 0.8 mm (10 Pack)
Ref: JC9	J Clasp - 0.9 mm (10 Pack)
Ref: JC10	J Clasp - 1.0 mm (10 Pack)

Silicone Mixing Cups

Suitable for mixing dental materials such as cement and acrylic (not suitable for silicone).

Available in 3 sizes



Stock Code	Description
Ref: P7220	Silicone Mixing Cup - Small ø 35 mm
Ref: P7230	Silicone Mixing Cup - Medium ø 48 mm
Ref: P7240	Silicone Mixing Cup - Large ø 70 mm

Plaster Bowls

Flexible multi purpose mixing bowls



Stock Code	Description
Ref: PBXL	Plaster Bowl Extra Large ø 150 mm
Ref: PBL	Plaster Bowl Large ø 140 mm
Ref: PBM	Plaster Bowl Medium ø 120 mm
Ref: PBS	Plaster Bowl Small ø 100 mm

Emery Cloth



Very firm bonding of the corundum to the base material makes this emery cloth very effective and durable.

Width 12 mm
Length 50 M



Stock Code	Description
Ref: D153012950	Emery Cloth (50 m roll)

Metroglue (Superglue)



Cyanoacrylate Adhesive

General purpose super glue.
Bonds in seconds. 20 g



Stock Code	Description
Ref: HGLUE	Metroglue Adhesive 20g

Metroblast - Alumina

Blasting Material

Ultra pure, Iron Free Aluminium Oxide for blast polishing with no surface metal loss. Metrodent blasting material is a white fused alumina manufactured by melting a special grade of calcined alumina (aluminium oxide) in an electric arc furnace using the block method; it is ultra pure aluminium oxide. An iron-free fine abrasive/refractory material for the toughest demands.

Advantages

- ✓ Ultra pure, Iron free aluminium oxide
- ✓ No surface metal loss



Stock Code	Description
Ref: MBAO250/10	Metroblast Alumina 250 μm - 10 Kg
Ref: MBAO110/10	Metroblast Alumina 110 μm - 10 Kg
Ref: MBAO50/10	Metroblast Alumina 50 μm - 10 Kg

Metroblast - Glass Beads

Blasting Material

In the manufacture of the glass beads alkali-lime glass is broken and prescreened. Subsequently, the glass fragments are formed in a furnace in the heat flow to balls. Glass beads are a mineral abrasive.

Advantages

- ✓ Glass beads for blast polishing
- ✓ No surface metal loss



Stock Code	Description
Ref: MBGB50/10	Metroblast Glass Beads 50 μm - 10 Kg

Steribim

Germicidal Pumice Compound

Destroys all germs and the germicidal effect lasts down to the last gram. Special additives make it kind to the skin and prevent the pumice from going sour in the polishing trough. Dust free.



Stock Code	Description
Ref: STB12	Steribim Germicidal Pumice 12 Kg

Metrotray

Light Curing Special Tray Material

For making customised individual trays for full / part posthetics or constructing base plates.



Stock Code	Description
Ref: DP08-00900	Pink - 50 Blanks
Ref: DP08-00901	Blue - 50 Blanks

Shellac Special Tray Blanks



Special Tray blanks are a firm favourite with Dentists and Technicians. Easily adapted to produce individual custom trays for the most accurate impressions.

Colour: Green
72 per box



Stock Code	Description
Ref: PAT112	Shellac Special Tray Blanks - Plain Upper (2.25 mm deep)
Ref: PAT122	Shellac Special Tray Blanks - Perforated Upper (2.25 mm deep)

Shellac Base Plates



Use Shellac Base plates when a very rigid base is needed to ensure maximum accuracy at the try-in stage of a denture.

Colour: Pink
100 per box



Stock Code	Description
Ref: BPL003	Shellac Base Plate - Upper (1.25 mm deep)
Ref: BPL004	Shellac Base Plate - Lower (1.25 mm deep)

Metrosil Plus



2 Part Silicone Duplicating Material

Metrosil is a two part addition cured silicone for the accurate duplication of master models. Mixed in a 1:1 ratio the material is blended to a uniform colour before pouring. Compared to reversible gels, silicone duplication has always been considered expensive. However, the growth of private and implant work has placed demands on the dental laboratory for high levels of accuracy and precision. Metrosil has excellent dimensional stability at an economical price.

1Kg Base
1Kg Catalyst



Stock Code	Description
Ref: MDSPLUS	Metrosil Plus (1Kg Set)



Metro Lab Putty



Silicone Modelling Material

Supplied as a base material and catalyst it is easily mixed by hand. Lab putty is routinely used for the production of copy denture moulds and denture repairs. Increasingly lab putty is being used for processing denture repairs. After sticking the parts of the broken denture together the repair model can be moulded quickly into the denture and placed on the bench to set. The fit surface is accurately reproduced before and after processing. The use of cold mould seal is not required. If for any reason the repair cannot be processed immediately, there is no danger of bubbles appearing in the cold cure. When repairing a midline crack it is advisable to make a complete model rather than one just over the crack as the denture can sometimes fracture when grinding out the repair. Lab putty is also used for matrices in the construction of partial dentures, for positioning backings, cold curing CoCr saddles etc.



Stock Code	Description
Ref: MLP	Metro Lab Putty Base - 5 Kg
Ref: MLPC	Metro Lab Putty Activator (Catalyst) - 45g

Laminated EVA Blanks

Tye Dye Range - 4mm


2 layer soft EVA laminated for strength.
Available in Ø 120 mm round and 125 mm square



Flamingo



Shamrock



Wilderness



Patriot



Midnight



Rainbow

Stock Code	Description
Ref: DR7954000R2	Patriot Round (12 Pack)
Ref: DR7924000R2	Wilderness Round (12 Pack)
Ref: DR7944000R2	Midnight Round (12 Pack)
Ref: DR7965000R2	Shamrock Round (12 Pack)
Ref: DR7934000R2	Flamingo Round (12 Pack)
Ref: DR7911111R2	Rainbow Round (12 Pack)
Ref: DR7954000	Patriot Square (12 Pack)
Ref: DR7924000	Wilderness Square (12 Pack)
Ref: DR7944000	Midnight Square (12 Pack)
Ref: DR7965000	Shamrock Square (12 Pack)
Ref: DR7934000	Flamingo Square (12 Pack)
Ref: DR7911111	Rainbow Square (12 Pack)

Laminated EVA Blanks

Glitter Range - 4mm


2 layer soft EVA laminated for strength.
Available in Ø 120 mm round and 125 mm square



Blue Glitter



Yellow Glitter



Purple Glitter



Pink Glitter



Orange Glitter



Green Glitter

Stock Code	Description
Ref: DR7965010R2	Blue Round (12 Pack)
Ref: DR7965030R2	Yellow Round (12 Pack)
Ref: DR7965050R2	Purple Round (12 Pack)
Ref: DR7965060R2	Pink Round (12 Pack)
Ref: DR7965020R2	Orange Round (12 Pack)
Ref: DR7965040	Green Round (12 Pack)
Ref: DR7965010	Blue Square (12 Pack)
Ref: DR7965030Y	Yellow Square (12 Pack)
Ref: DR7965050	Purple Square (12 Pack)
Ref: DR7965060	Pink Square (12 Pack)
Ref: DR7965020	Orange Square (12 Pack)
Ref: DR7965040	Green Square (12 Pack)

Laminated EVA Blanks

Fun Range - 4mm

2 layer soft EVA laminated for strength.
Available in Ø 120 mm round and 125 mm square



Blue Polkadot



Red Polkadot



White Polkadot



The Diagonal



Fangs




Dripping Blood

Stock Code	Description
Ref: DR9598079R2	Polka Dots - Blue Round (12 Pack)
Ref: DR9600929R2	Polka Dots - Red Round (12 Pack)
Ref: DR9601039R2	Polka Dots - White Round (12 Pack)
Ref: DR9600006R2	The Diagonal Round (12 Pack)
Ref: DR9599075R2	Fangs Round (12 Pack)
Ref: DR9599077R2	Dripping Blood Round (12 Pack)
Ref: DR9598079	Polka Dots - Blue Square (12 Pack)
Ref: DR9600929R	Polka Dots - Red Square (12 Pack)
Ref: DR9601039	Polka Dots - White Square (12 Pack)
Ref: DR9600006	The Diagonal Square (12 Pack)
Ref: DR9599075	Fangs Square (12 Pack)
Ref: DR9599077	Dripping Blood Square (12 Pack)


Laminated EVA Blanks

Classic Two Colour Range - 4mm


2 layer soft EVA laminated for strength.
Available in Ø 120 mm round and 125 mm square




Red/White




Purple/Gold



Blue/Gold




Green/Blue



Black/Red




Green/White




Green/Gold




White/Purple




Black/Gold



Blue/Green



Black/White




Blue/White

Stock Code

Description

- Ref: DR9598090R2 Black/Red Round (12 Pack)
- Ref: DR9598110R2 Black/Yellow (Gold) Round (12 Pack)
- Ref: DR9598130R2 Red/White Round (12 Pack)
- Ref: DR9598150R2 Blue/Yellow (Gold) Round (12 Pack)
- Ref: DR9598190R2 Green/Yellow (Gold) Round (12 Pack)
- Ref: DR9598210R2 Blue/Red Round (12 Pack)
- Ref: DR9598240R2 Black/White Round (12 Pack)
- Ref: DR9598260R2 Blue White Round (12 Pack)
- Ref: DR9598300R2 Maroon/Yellow (Gold) Round (12 Pack)
- Ref: DR9598340R2 Purple/Yellow (Gold) Round (12 Pack)
- Ref: DR9598360R2 Green/Blue Round (12 Pack)
- Ref: DR9598380R2 Green/White Round (12 Pack)
- Ref: DR9598450R2 White/Purple Round (12 Pack)



Maroon/Gold

- Ref: P110-BD01 Black/Red Square
- Ref: P110-BD02 Black/Yellow (Gold) Square (12 Pack)
- Ref: P110-BD03 Red/White Square (12 Pack)
- Ref: P110-BD04 Blue/Yellow (Gold) Square (12 Pack)
- Ref: P110-BD06 Green/Yellow (Gold) Square (12 Pack)
- Ref: P110-BD07 Blue/Red Square (12 Pack)
- Ref: P110-BD09 Black/White Square (12 Pack)
- Ref: P110-BD10 Blue White Square (12 Pack)
- Ref: P110-BD12 Maroon/Yellow (Gold) Square (12 Pack)
- Ref: P110-BD14 Purple/Yellow (Gold) Square (12 Pack)
- Ref: P110-BD15 Green/Blue Square (12 Pack)
- Ref: P110-BD16 Green/White Square (12 Pack)
- Ref: P110-BD20 White/Purple Square (12 Pack)

Laminated EVA Blanks

Classic Single Colour Range - 4mm

2 layer soft EVA laminated for strength.
Available in Ø 120 mm round and 125 mm square



Stock Code	Description
Ref: DR9598730R2	Black Round (12 Pack)
Ref: DR9600800R2	Green Round (12 Pack)
Ref: DR9600850R2	Orange Round (12 Pack)
Ref: DR9600870R2	Purple Round (12 Pack)
Ref: DR9600910R2	Red Round (12 Pack)
Ref: DR9598060R2	Blue Round (12 Pack)
Ref: DR9601020R2	White Round (12 Pack)
Ref: DR9601040R2	Yellow Round (12 Pack)
Ref: DR9596920R2	Clear Round (12 Pack)
Ref: P110-BKO	Black Square (12 Pack)
Ref: P110-G00	Green Square (12 Pack)
Ref: P110-ORA	Orange Square (12 Pack)
Ref: P110-P00	Purple Square (12 Pack)
Ref: P110-RED	Red Square (12 Pack)
Ref: P110-B00	Blue Square (12 Pack)
Ref: P110-W00	White Square (12 Pack)
Ref: P110-Y00	Yellow Square (12 Pack)
Ref: P110-000	Clear Square (12 Pack)
Ref: DR9598840R2	Tri Colour 4mm Red/White/Blue Round (12 Pack)
Ref: P110-BT03	Tri Colour 4mm Red/White/Blue Square (12 Pack)

Coping Material

Semi Clear

Easy to use and with stands grinding pressures well. A good choice for a patient that has minor grinding problems or a first time splint user.



Stock Code	Description
Ref: P210-001	Coping Material - 0.5 mm Square (12 Pack)
Ref: P210-002	Coping Material - 0.75 mm Square (12 Pack)
Ref: DR9613480R2	Coping Material - 0.75 mm Round (12 Pack)

Base Plate Material

Pure high impact styrene

Low melt temperature

Colour: Pink




Stock Code	Description
Ref: P231-001	Base Plate Material - 1.5 mm Square (12 Pack)
Ref: DR9616120R2	Base Plate Material - 1.5 mm Round (12 Pack)
Ref: P231-002	Base Plate Material - 2.0 mm Square (12 Pack)
Ref: DR9616130R2	Base Plate Material - 2.0 mm Round (12 Pack)
Ref: P231-003	Base Plate Material - 2.5 mm Square (12 Pack)

Splint Material

Manufacture of hard splints

This is a stiff clear material which is ideal for the manufacture of hard splints, thin splints, surgical trays and orthodontic retainers. The material will bond to acrylics to facilitate the manufacture of various dental appliances.



Stock Code	Description
Ref: P221-003	Splint Material 0.5 mm Square (12 Pack)
Ref: P221-006	Splint Material 1.0 mm Square (12 Pack)
Ref: P221-001	Splint Material 1.5 mm Square (12 Pack)
Ref: P221-002	Splint Material 2.0 mm Square (12 Pack)
Ref: P221-005	Splint Material 2.5 mm Square (12 Pack)
Ref: DR9615000R2	Splint Material 0.5 mm Round (12 Pack)
Ref: DR9640961R2	Splint Material 1.0 mm Round (12 Pack)
Ref: DR9640962R2	Splint Material 1.5 mm Round (12 Pack)
Ref: DR9640963R2	Splint Material 2.0 mm Round (12 Pack)
Ref: DR9640964R2	Splint Material 2.5 mm Round (12 Pack)

Crown and Bridge Material

Manufacture of temporary crowns

Thin clear rigid material 0.5 mm thick which is used for the manufacture of temporary bridges.



Stock Code	Description
Ref: DR9613680R2	Crown and Bridge Material - Round 0.5 mm (12 Pack)
Ref: P215-001	Crown and Bridge Material - Square 0.5 mm (12 Pack)

Dual Laminates

Two materials formed into a single laminated sheet. It was made especailly for use in bruxing cases. The soft liner assures patient comfort and the hard surface assures long-lasting wear.



Stock Code	Description
Ref: DR9603580R2	Dual Lamiates - Round Ø 120 x 3 mm (12 Pack)
Ref: P130-000	Dual Lamiates - Square 125 x 125 x 3 mm (12 Pack)

Nightguards

Made specifically for night time bruxers. The layering gives the thermoplastic a more abrasion-resistant characteristic that will withstand grinding better than traditional soft sheet materials.



Stock Code	Description
Ref: DR9602550R2	Nightguards - Round Ø 120 x 3 mm (12 Pack)
Ref: P120-000	Nightguards - Square 125 x 125 x 3 mm (12 Pack)

Wheel Polishers

Mounted HP Shank



Cotton Mop Mounted



Calico Wheel Mounted



Chamois 4 Fold Mounted



Stock Code	Description
Ref: HP/LW/CC	Cotton Mop Mounted - ø 21.0 mm
Ref: HP/LW/CAL	Calico Wheel Mounted - ø 20.0 mm
Ref: HP/LW/CH	Chamois 4 Fold Mounted - ø 21.0 mm

Lathe Brushes

Supplied in packs of 12



Black Calico Interleaved



Goat Hair Extra Soft



Stiff Bristle Nylon



Stiff Bristle

Stock Code	Description
Ref: 144C/P	Black Calico Interleaved (Plastic Hub) - ø 48.0 mm
Ref: 144G/P	Goat Hair Extra Soft (Plastic Hub) - ø 48.0 mm
Ref: 36A/NYLON	Stiff Bristle Nylon (Plastic Hub) - ø 48.0 mm
Ref: 36A	Stiff Bristle (Plastic Hub) - ø 48.0 mm

Mounted Cups



HP and RA Shanks



HP Stiff Bristle -
White



HP Stiff Nylon -
White



RA Stiff Bristle -
White



RA Stiff Nylon -
White

Stock Code	Description
Ref: C/WB10HP	HP Stiff Bristle White Mounted Cup (10 Pack) - ø 10.0 mm
Ref: C/WN10HP	HP Stiff Nylon White Mounted Cup (10 Pack) - ø 10.0 mm
Ref: JC/WBRA	RA Stiff Bristle White Mounted Cup (100 Pack) - ø 10.0 mm
Ref: JC/WNRA	RA Stiff Nylon White Mounted Cup (100 Pack) - ø 10.0 mm

Rubber Prophy Polishing Cups



Supplied in packs of 100



Stock Code	Description
Ref: TPY002	RA 4 Web Latch Rubber Prophy Polishing Cup (100 Pack)
Ref: TPR004A	RA 4 Web Snap On Rubber Prophy Polishing Cup (100 Pack)

Polishing Mops

Supplied in Packs 10



Calico Mop -
25 Fold Stiff



Calico Mop -
25 Fold Medium



Calico Mop -
30 Fold Soft Weave




Wool Mop -
For Lathe

Stock Code

Description

Ref: H/4/25
Ref: M/4/25
Ref: F/4/30
Ref: BR60L

Calico Mop 25 Fold Stiff (10 Pack)- ø 100.0 mm
Calico Mop 25 Fold Medium(10 Pack) - ø 100.0 mm
Calico Mop 30 Fold Soft Weave(10 Pack) - ø 100.0 mm
Wool Mop For Lathe(10 Pack) - ø 100.0 mm



Lathe Brushes

Supplied in Packs 10



Calico Interleaved
Lathe Brush



Straight
Lathe Brush



Goat Hair Extra Soft
Lathe Brush

Stock Code	Description
Ref: 86P	Calico Interleaved Lathe Brush (10 Pack)- ø 80.0 mm
Ref: 8P	Straight Lathe Brush (10 Pack)- ø 80.0 mm
Ref: 150	Goat Hair Extra Soft Lathe Brush (10 Pack)- ø 80.0 mm


Lathe Brushes

Supplied in Packs 10



Black Extra Stiff - 3 Row
Lathe Brush

Black Extra Stiff - 2 Row
Lathe Brush



Stock Code	Description
Ref: 22P	Black Extra Stiff 3 Row Lathe Brush (10 Pack)- ø 63.0 mm
Ref: 29A	Black Extra Stiff 2 Row Lathe Brush (10 Pack)- ø 50.0 mm

Spatulas for Plaster, Alginate and Elastomers



Stock Code	Description
Ref: PKL	1. Plaster Knife - Large
Ref: PKM	2. Plaster Knife - Medium with flask opener 1
Ref: PKS	3. Plaster Knife - Small
Ref: SPT2	4. Spatula - Large flexible
Ref: SPT3	5. Spatula - Medium flexible
Ref: SPT4	6. Spatula - Small flexible
Ref: SPT5	7. Spatula - Small rigid
Ref: SPT6	8. Spatula for plaster - Large rigid

Wax Knives



Stock Code	Description
Ref: WKS	1. Fahnenstock - Small Wax Knife
Ref: WKS	2. Fahnenstock - Large Wax Knife



Le Cron Carver

Modelling Instrument

Wax and modelling carver



Stock Code	Description
Ref: WMD	Le Cron Wax and Modelling Carver



No.5 Wax Carver

Modelling Instrument

Length 170 mm



Stock Code	Description
Ref: WMC	No.5 Wax Spatula

Calcium Carbonate, Percipitated



Whiting Powder

HS Tariff: 2836.50.00



Stock Code	Description
Ref: RMS/250	Calcium Carbonate, Percipitated (Whiting Powder) - 25 Kg

Prophylaxis Paste

Original

200g Jar



Stock Code	Description
Ref: LAB030	Prophylaxis Paste - 200g (12 Jars)

Flange Mouthpacks



Non Sterile

Absorbent sponges to pack the airway in order to prevent inhalation/ ingestion of debris, etc.

Autoclave before use.



Stock Code	Description
Ref: MPACK	Flange Mouth Pack (12 Pack)

Bur Cleaner



With Sliding Brush



Stock Code	Description
Ref: BUR/002	Bur Cleaning Brush With Slide (12 Pack)

Shipping Instructions

In order to provide you with a prompt and most effective service possible for seafreight or airfreight consignments we recommend that your purchasing orders should allow for partial shipments and indicate the following details:

1. Method of shipment required
2. Delivery location
3. Port/airport of discharge
4. Shipping marks required
5. Shipping documents required
6. Your banking details

When shipping method is not specified, the mode of shipment used will be based on the nature of goods and the size of the consignment. Airfreight consignments will be made upon specific request with the exception of hazardous goods.

Handling of Hazardous Goods

Before any consignment of hazardous goods is made, we ensure that proper classification, packing and labelling are carried out in accordance with the international regulations governing the handling and transportation of hazardous goods. Packing and handling charges may be added to defray the cost incurred. Metrodent believes the information contained in this catalogue is accurate. However the information is subject to change without notice and should not be construed as a commitment by Metrodent Limited. Metrodent assumes no responsibility for any errors that may appear in this catalogue.

Loss and Damage in Transit

Upon receipt of the consignment, we advise you strongly to check for damages and/or shortages. In the event that shortages and/or damages have been ascertained, you should without delay distraint the carrier to draw-up a report of findings. This document is a mandatory part of the claim to be made on the carrier. You should also notify the carrier within 14 days from date of delivery of any concealed damages and/or shortages.

Return of Goods for Credit

To return goods for credit you must first obtain a written permission and shipping instructions from Metrodent Limited. New and unused goods packed in their original packages may be returned for credit within 90 days of date of shipment. Metrodent Limited reserve the rights to charge a restocking fee equal to up to 25% of the selling price. Authorised goods returned for credit should be shipped insured and freight prepaid.

Prices

Our prices are deemed to be ex-works, exclusive of packing and handling costs and all other related charges (freight, insurance, documentation etc.)



The information and prices are subject to change without notice and should not be construed as a commitment by Metrodent Limited.
Metrodent assumes no responsibility for any errors that may appear in this product guide.